

## HIGH DEFINITION Still Capture

Capture, print and export 1080p still images with the HD Still Capture System from Conmed Linvatec. This HD Still Capture System is simple, easy-to-use and provides a range of image-preview options for the discerning surgeon. With its LINUX operating system and embedded system design, the Still Capture System starts up fast, recovers elegantly from power failures and is more immune to computer viruses and network threats.

- TRUE-HD 1080p or 720p still capture
- Full screen or thumbnail display of captured HD images
- Full image review and delete capability
- Print, save to USB or send images to PACS via DICOM
- Easy and intuitive user interface
- Upgradeable FPGA platform



HD Image  
Capture



Print



DICOM



USB Flash  
Memory




\*HD Monitor  
sold separately



# HD 1080p

# HIGH DEFINITION Still Capture



## HD 1080p

### VP1600 HD Still Capture System



1-2-3 simple and easy to use

### Related ConMed Linvatec Products



HD 1080p  
Monitor

IM4000  
HD Camera

#### Catalog Number

VP1600

#### Description

HD Still Capture System

Please refer to the VP1600 Operator's Manual for a full list of compatible accessories such as printer, USB storage, cables etc.

#### Specifications

Operating System:	LINUX
Video Input/Output:	DVI-D, HD Video (1080p or 720p)
Image Matrix	
1080p:	1920 x 1080
720p:	1280 x 720

#### Control

Still Capture	3.5mm Mini Phono x 2
Nurses Assistant®:	RS-232
Data I/O	
Keyboard:	USB 2.0 Type 1
Printer:	USB 2.0 Type 1
USB Storage Device	USB 2.0 Type 1 (Front or Rear)

#### Typical Operating Requirements

Input Voltage:	100 – 240 VAC
Input Line Frequency:	50/60 Hz
Power Consumption:	100 VA maximum
Fuse:	(2X) T 1.0 A, 250 V

#### Dimensions

Controller:	3.63" (H) x 14.0" (W) x 14.0" (L) 9.2cm (H) x 35.6cm (W) x 35.6cm (L)
Keyboard:	0.9" (H) x 11.1" (W) x 5.3" (D) 2.3cm (H) x 28.2cm (W) x 13.5cm (D)
Weight:	7 lbs. (3.2 kg)

#### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard  
Largo, FL 33773-4908  
(727) 392-6464  
Customer Service: 1-800-237-0169  
FAX: (727) 399-5256  
International FAX: +1 (727) 397-4540  
email: customer\_service@linvatec.com  
www.linvatec.com