


Large and Small Bone Catalog

 **smith&nephew**
DYONICS[®] POWER
Battery Instruments and Disposables

When performance is everything,
trust the power of DYONICS.





Battery Instruments

3 **Battery Instruments**

4 Large Bone Battery Instruments & Accessories

6 Small Bone Battery Instruments & Accessories

Disposables

8 **Overview of Blade Styles**

9 **Large Bone Blades**

10 Oscillating Saw Blades

19 Reciprocating Saw Blades

20 Knee Implant Cross Reference

21 **Small Bone Blades**

22 ACL Saw Blades

22 Sagittal Saw Blades

31 Reciprocating Saw Blades

33 **Rasps, Twist Drills, K-Wires, Steinmann Pins & Pin Protectors**

34 Rasps

35 Twist Drills

37 K-Wires

38 Steinmann Pins

39 Pin Protectors





Battery Instruments

Large Bone

- 4 Large Bone Battery Instruments
- 5 Accessories

Small Bone

- 6 Small Bone Battery Instruments
- 7 Accessories

Cordless Instruments

Depend on DYONICS® POWER to deliver the performance you demand with reliability and quality you can count on.



Oscillating Saw

- Dual settings for enhanced speed control at 12,000 and 15,000 cpm
- Rotatable cutting head aids cutting at various angles
- Secure and easy to operate latch-style locking system



Drill / Reamer

- Commanding 750 rpm drilling power with 50 inch-pound of torque
- Effective 250 rpm reaming speed with 150 inch-pound of torque



Reciprocating Saw

- Fast-moving blades at 15,000 cpm can be rotated a full 360°
- Twist-lock blade hub provides efficient and easy removal or exchange



Precision molded from tough, lightweight machined aluminum

Ergonomically efficient handpiece distributes and balances weight to help reduce muscle tension while in the hand

Reliable and long-lasting NiMH rechargeable battery simplifies battery management

Oscillating Saw

DYONICS POWER Large Bone Battery Instrumentation

| Reference # | Description |
|-------------|--|
| 72204091 | DYONICS POWER Large Bone Drill/Reamer |
| 72204092 | DYONICS POWER Large Bone Oscillating Saw |
| 72204093 | DYONICS POWER Large Bone Reciprocating Saw |

Power Supplies

Powerful, 14.4 volt Nickel Metal Hydride (NiMH) technology batteries offer long run time and exceptional charge-holding ability.

- Sterilization options include autoclavable or aseptic
- No heavy metals means no special disposal requirements
- Full one-year warranty* can help minimize ownership cost



Universal Large Battery



Universal Battery Charger (four port)



Aseptic Battery

DYONICS POWER Batteries and Charging Stations

| Reference # | Description |
|-------------|--|
| 72204113* | Large NiMH Battery Pack (autoclavable) |
| 72204109 | Universal Battery Charger (four port) |
| 72204105 | Aseptic Battery Housing |
| 72204106 | Aseptic NiMH Battery Pack |
| 72204107 | Aseptic Battery Charging Adapter |
| 72204108 | Aseptic Battery Transfer Shield |

*One-year warranty applies to 72204113 Large NiMH Battery Pack (autoclavable)

Accessories

Large Bone Drilling and Reaming Couplers

Drilling and reaming couplers offer exceptional speed and torque to tackle a wide variety of surgical procedures.



Drilling Couplers

| Reference # | Description |
|-------------|--|
| 72204095 | Jacobs-style Coupler 1/4" with Key |
| 72204096 | Keyless Coupler 1/4" |
| 72204097 | Keyless Coupler 1/8" |
| 72204098 | Synthes Quick-connect Coupler |
| 72204102 | Zimmer (Hudson Modified Trinkle) Coupler |
| 72204103 | Hudson Coupler |



Reaming Couplers

| Reference # | Description |
|-------------|---|
| 72204099 | Hudson High-torque Reamer Coupler |
| 72204094 | Zimmer (Hudson Modified Trinkle) Reamer Coupler |
| 72204100 | Jacobs-style High-torque Reamer Coupler 1/4" |
| 72204101 | Harris/Aesculap Reamer Coupler |
| 72204104 | Synthes Reaming Coupler |

Large Bone Sawing

Available with quick-connect keyless connections or locking keys, sagittal and reciprocating saws provide precision and tough cutting.




Sawing

| Reference # | Description |
|-------------|--|
| 72204117 | Sagittal Saw Coupler with Key |
| 72204118 | Sagittal Saw Coupler Keyless |
| 72204119 | Heavy-duty Sagittal Saw Coupler with Key |
| 72204115 | Reciprocating Saw Coupler |

Large Bone Wire & Pin Driving

Automatic wire and pin collectors facilitate ease of placement, and when used with the oscillate feature, accurate placement can be achieved without spooling of soft tissue.




Wire & Pin Driving

| Reference # | Description |
|-------------|-----------------|
| 72204111 | Wire/Pin Driver |
| 72204112 | Pin Driver |

Large Bone Sterilization

| Reference # | Description |
|-------------|--------------------|
| 72204110 | Sterilization Case |



Dual-trigger Cordless Handpiece

Depend on DYONICS® POWER to deliver the performance you demand with reliability and quality you can count on.



DYONICS POWER Small Bone Battery Instrumentation

| Reference # | Description |
|-------------|---------------------------------------|
| 72204114 | DYONICS POWER Small Battery Handpiece |

BATTERY INSTRUMENTS

Power Supplies

Powerful, 14.4 volt Nickel Metal Hydride (NiMH) technology batteries offer long run time and exceptional charge-holding ability.

- Flexible battery options are available in two sizes
- Sterilization options include autoclavable or aseptic
- No heavy metals means no special disposal requirements
- Full one-year warranty* can help minimize ownership cost

DYONICS POWER Batteries and Charging Stations

| Reference # | Description |
|-------------|--|
| 72204125* | Small NiMH Battery Pack (autoclavable) |
| 72204113* | Large NiMH Battery Pack (autoclavable) |
| 72204109 | Universal Battery Charger (four port) |
| 72204105 | Aseptic Battery Housing |
| 72204106 | Aseptic NiMH Battery Pack |
| 72204107 | Aseptic Battery Charging Adapter |
| 72204108 | Aseptic Battery Transfer Shield |

*One-year warranty applies to 72204125 Small NiMH Battery Pack (autoclavable) and 72204113 Large NiMH Battery Pack (autoclavable)



Small Battery



Universal Large Battery



Aseptic Battery



Universal Battery Charger (four port)

Accessories

Small Bone Drilling & Reaming Couplers

Drilling and reaming couplers offer the speed and torque necessary to tackle a wide variety of surgical procedures.



72204122



72204120



72204124



72204602

Drilling Couplers

| Reference # | Description |
|-------------|-----------------------------|
| 72204122 | 1/4" Jacobs Drill |
| 72204120 | 5/32" Jacobs Drill |
| 72204124 | Synthes Quick Connect Drill |
| 72204602 | 1/4" Keyless Jacobs Coupler |



72204121



72204116



72204123

Reaming Couplers

| Reference # | Description |
|-------------|---|
| 72204121 | 1/4" Jacobs High Torque Reamer |
| 72204116 | Hudson High Torque Reamer |
| 72204123 | Zimmer Reamer (Hudson Modified Trinkle) |

Small Bone Sterilization

Sterilization Case

| Reference # | Description |
|-------------|-----------------------------------|
| 72204126 | Sterilization Case, Small Battery |



Small Bone Wire & Pin Driving

Automatic wire and pin collectors facilitate ease of placement, and when used with the oscillate feature, accurate placement can be achieved without spooling of soft tissue.



72204127



72204128

Wire & Pin Drivers

| Reference # | Description |
|-------------|--|
| 72204127 | Wire Driver – accepts wires from .028" (0.7 mm)–.062" (1.6 mm) |
| 72204128 | Pin Driver – accepts pins from .078" (1.7 mm)–.125" (3.2 mm) |

Small Bone Sawing

Available with quick-connect keyless connections or locking keys, sagittal and reciprocating saws provide precision and tough cutting.



72204117



72204118



72204119



72204115

Sawing

| Reference # | Description |
|-------------|----------------------------------|
| 72204117 | Sagittal Saw with Key |
| 72204118 | Keyless Sagittal Saw |
| 72204119 | Heavy-duty Sagittal Saw with Key |
| 72204115 | Reciprocating Saw |


Cut Face Reference

Ideal for a variety of arthroscopic procedures, unique blade designs are available to suite the specific needs of the surgeon and the clinical demands of any procedure.




StableCut®

- Stable and straight cutting
- Patented tooth design
- Efficient debris ejection




Diamond Ground

- Angles are ground into teeth
- Two rows of teeth
- Cutting thickness same as blade



Dual-Cut

- Small double teeth
- Aligned on straight plane
- Channels help move debris from cutting surface




Aggressive Offset

- Aggressive cutting
- Teeth are bent out forming two rows
- Comparable thickness to PhotoEtch Aggressive



PhotoEtch Aggressive

- Aggressive cutting
- Photo-chemically etched dual thickness
- Reduced friction when used with a slotted cutting guide



Standard Offset

- Teeth are alternately bent outwards
- Controlled cutting
- Geometry similar to Diamond Ground



Large Bone Blades

Oscillating Saw Blades

- 10 StableCut
- 11 Dual-Cut
- 12 PhotoEtch Aggressive
- 14 Diamond Ground
- 15 Aggressive Offset
- 17 Standard Offset

Reciprocating Saw Blades


- 19 Reciprocating Saw Blades

Large Bone Blades




Oscillating Saw Blades - StableCut®

▶ Blades are sterile and packaged 6 per box




Part Number: **72204192**

W 13 mm CT 1.0 mm TPI 12
UL 90 mm MT 1.0 mm




Part Number: **72204188**

W 13 mm CT 1.27 mm TPI 12
UL 90 mm MT 1.27 mm



Part Number: **72204164**

W 19 mm CT 1.35 mm TPI 12
UL 90 mm MT 1.35 mm




Part Number: **72204178**

W 21.3 mm CT 1.27 mm TPI 12
UL 90 mm MT 1.27 mm



Part Number: **72204189**

W 21.3 mm CT 1.19 mm TPI 12
UL 90 mm MT 1.19 mm



Part Number: **72204194**

W 25 mm CT 10 mm TPI 12
UL 90 mm MT 10 mm


DISPOSABLES

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |




Oscillating Saw Blades - StableCut®

▶ Blades are sterile and packaged 6 per box




Part Number: **72204175**

| | | |
|-----------------|-------------------|---------------|
| W 25 mm | CT 1.27 mm | TPI 12 |
| UL 90 mm | MT 1.27 mm | |




Part Number: **72204170**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.47 mm | TPI 12 |
| UL 90 mm | MT 1.47 mm | |




Oscillating Saw Blades - Dual-Cut

▶ Blades are sterile and packaged 6 per box




Part Number: **72204193**

| | | |
|-----------------|-------------------|---------------|
| W 18 mm | CT 0.89 mm | TPI 16 |
| UL 90 mm | MT 0.89 mm | |



Part Number: **72204186**

| | | |
|-----------------|-------------------|---------------|
| W 18 mm | CT 1.27 mm | TPI 16 |
| UL 90 mm | MT 1.27 mm | |




Part Number: **72204191**

| | | |
|-----------------|-------------------|---------------|
| W 18 mm | CT 1.35 mm | TPI 16 |
| UL 90 mm | MT 1.35 mm | |




Oscillating Saw Blades - Dual-Cut

▶ Blades are sterile and packaged 6 per box




Part Number: **72204190**

W 25 mm CT 1.27 mm TPI 16
 UL 90 mm MT 1.27 mm




Part Number: **72204176**

W 25 mm CT 1.27 mm TPI 16
 UL 100 mm MT 1.27 mm




Part Number: **72204173**

W 25 mm CT 1.35 mm TPI 16
 UL 90 mm MT 1.35 mm




Part Number: **72204185**

W 25 mm CT 1.47 mm TPI 16
 UL 90 mm MT 1.47 mm



Oscillating Saw Blades - PhotoEtch Aggressive


▶ Blades are sterile and packaged 6 per box



Part Number: **72204145**


W 19.5 mm CT 1.27 mm TPI 10
 UL 90 mm MT 1.27 mm

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |




Oscillating Saw Blades - PhotoEtch Aggressive

▶ Blades are sterile and packaged 6 per box




Part Number: **72204151**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 10 |
| UL 90 mm | MT 0.89 mm | |




Part Number: **72204154**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 1.20 mm | TPI 10 |
| UL 90 mm | MT 1.20 mm | |




Part Number: **72204147**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 1.27 mm | TPI 10 |
| UL 95 mm | MT 1.27 mm | |




Part Number: **72204163**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 0.90 mm | TPI 10 |
| UL 90 mm | MT 0.89 mm | |



Part Number: **72204148**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.27 mm | TPI 10 |
| UL 90 mm | MT 1.27 mm | |



Part Number: **72204182**


| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.37 mm | TPI 10 |
| UL 90 mm | MT 1.37 mm | |

DISPOSABLES




Oscillating Saw Blades - Diamond Ground

▶ Blades are sterile and packaged 6 per box




Part Number: **72204180**

| | | |
|-----------------|-------------------|---------------|
| W 10 mm | CT 1.27 mm | TPI 12 |
| UL 85 mm | MT 1.27 mm | |




Part Number: **72204152**

| | | |
|-----------------|-------------------|---------------|
| W 13 mm | CT 1.19 mm | TPI 12 |
| UL 90 mm | MT 1.19 mm | |




Part Number: **72204166**

| | | |
|------------------|-------------------|---------------|
| W 19.0 mm | CT 1.35 mm | TPI 12 |
| UL 85 mm | MT 1.35 mm | |




Part Number: **72204183**

| | | |
|------------------|-------------------|---------------|
| W 19.1 mm | CT 1.27 mm | TPI 12 |
| UL 85 mm | MT 1.27 mm | |



Part Number: **72204156**


| | | |
|------------------|-------------------|---------------|
| W 19.1 mm | CT 1.19 mm | TPI 12 |
| UL 105 mm | MT 1.19 mm | |



Part Number: **72204171**


| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.19 mm | TPI 12 |
| UL 90 mm | MT 1.19 mm | |

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |




Oscillating Saw Blades - Diamond Ground

▶ Blades are sterile and packaged 6 per box




Part Number: **72204146**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.20 mm | TPI 12 |
| UL 90 mm | MT 1.20 mm | |




Oscillating Saw Blades - Aggressive Offset

▶ Blades are sterile and packaged 6 per box




Part Number: **72204150**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 18 |
| UL 71 mm | MT 0.76 mm | |




Part Number: **72204158**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 10 |
| UL 71 mm | MT 0.76 mm | |



Part Number: **72204177**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 1.40 mm | TPI 18 |
| UL 71 mm | MT 1.19 mm | |



Part Number: **72204181**


| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 10 |
| UL 80 mm | MT 0.76 mm | |

DISPOSABLES




Oscillating Saw Blades - Aggressive Offset

▶ Blades are sterile and packaged 6 per box




Part Number: **72204160**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 1.40 mm | TPI 10 |
| UL 80 mm | MT 1.19 mm | |




Part Number: **72204159**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 18 |
| UL 86 mm | MT 0.76 mm | |




Part Number: **72204187**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 0.90 mm | TPI 10 |
| UL 90 mm | MT 0.76 mm | |




Part Number: **72204184**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.40 mm | TPI 10 |
| UL 90 mm | MT 1.19 mm | |



Part Number: **72204179**

| | | |
|------------------|-------------------|---------------|
| W 33.3 mm | CT 0.90 mm | TPI 18 |
| UL 63 mm | MT 0.64 mm | |



Part Number: **72204167**


| | | |
|------------------|-------------------|---------------|
| W 33.5 mm | CT 0.90 mm | TPI 18 |
| UL 64 mm | MT 0.76 mm | |

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |




Oscillating Saw Blades - Standard Offset

▶ Blades are sterile and packaged 6 per box



Part Number: **72204153**

| | | |
|-----------------|-------------------|---------------|
| W 11 mm | CT 1.10 mm | TPI 18 |
| UL 30 mm | MT 0.64 mm | |



Part Number: **72204169**

| | | |
|-----------------|-------------------|---------------|
| W 13 mm | CT 0.90 mm | TPI 18 |
| UL 64 mm | MT 0.64 mm | |




Part Number: **72204172**

| | | |
|-----------------|-------------------|---------------|
| W 13 mm | CT 1.50 mm | TPI 12 |
| UL 70 mm | MT 1.19 mm | |




Part Number: **72204165**

| | | |
|-----------------|-------------------|---------------|
| W 13 mm | CT 1.00 mm | TPI 16 |
| UL 90 mm | MT 0.89 mm | |



Part Number: **72204144**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.90 mm | TPI 18 |
| UL 63 mm | MT 0.64 mm | |




Part Number: **72204157**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.00 mm | TPI 12 |
| UL 90 mm | MT 0.89 mm | |



Oscillating Saw Blades - Standard Offset

▶ Blades are sterile and packaged 6 per box



Part Number: **72204174**

| | | |
|------------------|-------------------|---------------|
| W 25.4 mm | CT 1.40 mm | TPI 12 |
| UL 90 mm | MT 1.27 mm | |

| W | UL | CT | MT | TPI |
|-------|---------------|---------------|--------------------|----------------|
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |


Reciprocating Saw Blades

▶ Blades are sterile and packaged 6 per box




Part Number: **72204161**

| | | | | | |
|----|---------|----|---------|-----|----|
| W | 9.4 mm | CT | 1.00 mm | TPI | 18 |
| UL | 52.8 mm | MT | 0.64 mm | | |




Part Number: **72204168**

| | | | | | |
|----|---------|----|---------|-----|----|
| W | 12.2 mm | CT | 1.80 mm | TPI | 18 |
| UL | 80 mm | MT | 1.80 mm | | |




Part Number: **72204149**

| | | | | | |
|----|---------|----|---------|-----|----|
| W | 12.2 mm | CT | 1.00 mm | TPI | 18 |
| UL | 87.9 mm | MT | 0.64 mm | | |



Part Number: **72204155**

| | | | | | |
|----|---------|----|---------|-----|----|
| W | 12.7 mm | CT | 0.90 mm | TPI | 20 |
| UL | 73.5 mm | MT | 0.81 mm | | |



Part Number: **72204162**

| | | | | | |
|----|---------|----|---------|-----|----|
| W | 12.7 mm | CT | 0.90 mm | TPI | 16 |
| UL | 76 mm | MT | 0.90 mm | | |

DISPOSABLES

Knee Implant Cross Reference

▶ Blades are sterile and packaged 6 per box

| Implant | StableCut Blade | Dual-Cut Blade | Standard Blade | Thickness |
|---|--|--------------------|--|-----------|
| Aesculap Columbus® | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| Biomet AGC®, Ascent™, Maxim®, Vanguard® | | 72204173 – Page 12 | 72204182 – Page 13 | 1.37 mm |
| Biomet Performance® | 72204194 – Page 10 72204192 – Page 10 | | 72204151 – Page 13 72204157 – Page 17 | 1.07 mm |
| Corin All Knee Systems | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| DePuy AMK®, RHK, S-ROM® | | | 72204157 – Page 17 | 0.89 mm |
| DePuy Complete™, L.C.S.® | 72204170 – Page 11 | 72204185 – Page 12 | | 1.47 mm |
| DePuy P.F.C. Sigma® | | | 72204154 – Page 13 72204152 – Page 14 | 1.19 mm |
| DJO 3D-Knee™, Foundation® | 72204194 – Page 10 72204192 – Page 10 | | 72204151 – Page 13 72204157 – Page 17 | 1.00 mm |
| Exactech Optetrak® | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| Smith & Nephew TC-Plus°, RT-Plus° | 72204194 – Page 10 72204192 – Page 10 | 72204193 – Page 11 | 72204151 – Page 13 72204157 – Page 17 | 1.00 mm |
| Smith & Nephew GENESIS° | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| Smith & Nephew JOURNEY°, Profix° | 72204164 – Page 10 | 72204191 – Page 11 | 72204166 – Page 14 | 1.37 mm |
| Stryker All Knee Systems | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| Wright Advance®, Advantim® | 72204175 – Page 11 72204188 – Page 10 | 72204190 – Page 12 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |
| Wright Axiom® | | | 72204157 – Page 17 | 0.89 mm |
| Zimmer IB II, MG II® | | | 72204157 – Page 17 | 0.89 mm |
| Zimmer Natural Knee® | 72204194 – Page 10 72204192 – Page 10 | 72204193 – Page 11 | 72204151 – Page 13 72204157 – Page 17 | 1.00 mm |
| Zimmer Nexgen® | 72204175 – Page 11 72204188 – Page 10 | 72204186 – Page 11 | 72204145 – Page 12 72204183 – Page 14 | 1.27 mm |



Small Bone Blades

ACL Saw Blades

22 ACL Saw Blades

Sagittal Saw Blades


22 Sagittal Saw Blades

Reciprocating Saw Blades

31 Reciprocating Saw Blades

Sagittal Saw Blades - ACL

▶ Blades are sterile and packaged 6 per box




Part Number: **72204273**

W 8 mm **CT** 0.60 mm **TPI** 18

UL 12 mm **MT** 0.38 mm

Titanium nitride coated



Part Number: **72204277**


W 9.5 mm **CT** 0.60 mm **TPI** 18

UL 10 mm **MT** 0.38 mm

Titanium nitride coated

Sagittal Saw Blades - Keyless


▶ Blades are sterile and packaged 6 per box



Part Number: **72204236**

W 4.5 mm **CT** 0.6 mm **TPI** 36


UL 25.5 mm **MT** 0.4 mm



Part Number: **72204278**

W 5.5 mm **CT** 0.4 mm **TPI** 32

UL 9.5 mm **MT** 0.4 mm



Part Number: **72204247**


W 5.5 mm **CT** 0.4 mm **TPI** 32

UL 18 mm **MT** 0.4 mm

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |


Sagittal Saw Blades - Keyless

▶ Blades are sterile and packaged 6 per box




Part Number: **72204266**

| | | |
|-----------------|------------------|---------------|
| W 5.5 mm | CT 0.4 mm | TPI 32 |
| UL 22 mm | MT 0.4 mm | |




Part Number: **72204255**

| | | |
|-----------------|------------------|---------------|
| W 6 mm | CT 0.5 mm | TPI 34 |
| UL 34 mm | MT 0.4 mm | |




Part Number: **72204237**

| | | |
|-------------------|------------------|---------------|
| W 5.6 mm | CT 0.6 mm | TPI 18 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204283**

| | | |
|-----------------|------------------|---------------|
| W 7 mm | CT 0.4 mm | TPI 32 |
| UL 13 mm | MT 0.4 mm | |



Part Number: **72204257**

| | | |
|-----------------|------------------|---------------|
| W 7 mm | CT 0.4 mm | TPI 32 |
| UL 16 mm | MT 0.4 mm | |



Part Number: **72204254**


| | | |
|-------------------|------------------|---------------|
| W 7 mm | CT 0.4 mm | TPI 32 |
| UL 29.5 mm | MT 0.4 mm | |

DISPOSABLES

Sagittal Saw Blades - Keyless


▶ Blades are sterile and packaged 6 per box

DISPOSABLES




Part Number: **72204264**

W 9 mm CT 0.4 mm TPI 32
 UL 17 mm MT 0.4 mm




Part Number: **72204242**

W 9 mm CT 0.4 mm TPI 22
 UL 31 mm MT 0.4 mm




Part Number: **72204252**

W 9 mm CT 0.4 mm TPI 32
 UL 31 mm MT 0.4 mm




Part Number: **72204282**

W 9.4 mm CT 0.6 mm TPI 36
 UL 15 mm MT 0.4 mm



Part Number: **72204227**

W 9.5 mm CT 0.5 mm TPI 34
 UL 25.4 mm MT 0.4 mm




Part Number: **72204258**

W 9.5 mm CT 0.5 mm TPI 31
 UL 24.4 mm MT 0.4 mm

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |


Sagittal Saw Blades - Keyless

▶ Blades are sterile and packaged 6 per box




Part Number: **72204230**

| | | |
|-------------------|------------------|---------------|
| W 9.5 mm | CT 0.6 mm | TPI 18 |
| UL 25.4 mm | MT 0.4 mm | |



Part Number: **72204275**

| | | |
|-------------------|------------------|---------------|
| W 9.5 mm | CT 1.0 mm | TPI 34 |
| UL 25.4 mm | MT 0.8 mm | |




Part Number: **72204262**

| | | |
|-----------------|------------------|---------------|
| W 9.5 mm | CT 0.5 mm | TPI 34 |
| UL 41 mm | MT 0.4 mm | |



Part Number: **72204263**

| | | |
|-----------------|------------------|---------------|
| W 10 mm | CT 0.5 mm | TPI 34 |
| UL 35 mm | MT 0.6 mm | |



Part Number: **72204153**

| | | |
|-----------------|-------------------|---------------|
| W 11 mm | CT 1.1 mm | TPI 18 |
| UL 30 mm | MT 0.64 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).



Part Number: **72204172**


| | | |
|-----------------|-------------------|---------------|
| W 13 mm | CT 1.5 mm | TPI 12 |
| UL 70 mm | MT 1.19 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).

DISPOSABLES

Sagittal Saw Blades - Keyless


▶ Blades are sterile and packaged 6 per box



Part Number: **72204169**


| | | |
|-----------------|------------------|---------------|
| W 13 mm | CT 0.9 mm | TPI 18 |
| UL 64 mm | MT 0.9 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).



Part Number: **72204286**

| | | |
|-----------------|-------------------|---------------|
| W 14 mm | CT 1.00 mm | TPI 18 |
| UL 16 mm | MT 0.6 mm | 90° |




Part Number: **72204243**

| | | |
|-----------------|------------------|---------------|
| W 14 mm | CT 0.6 mm | TPI 18 |
| UL 41 mm | MT 0.4 mm | |



Part Number: **72204260**


| | | |
|------------------|------------------|---------------|
| W 19.5 mm | CT 0.6 mm | TPI 18 |
| UL 41 mm | MT 0.4 mm | |



Part Number: **72204144**

| | | |
|------------------|------------------|---------------|
| W 19.5 mm | CT 0.9 mm | TPI 18 |
| UL 63 mm | MT 0.6 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).



Part Number: **72204150**

| | | |
|------------------|------------------|---------------|
| W 19.5 mm | CT 0.9 mm | TPI 18 |
| UL 71 mm | MT 0.8 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |

Sagittal Saw Blades - Keyless


▶ Blades are sterile and packaged 6 per box



Part Number: **72204158**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 0.9 mm | TPI 10 |
| UL 71 mm | MT 0.76 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).




Part Number: **72204177**

| | | |
|------------------|-------------------|---------------|
| W 19.5 mm | CT 1.4 mm | TPI 18 |
| UL 71 mm | MT 1.19 mm | |

Note: Requires Heavy-duty Sagittal Saw Coupler with Key (REF 72204119, page 5).


Sagittal Saw Blades - Keyed

▶ Blades are sterile and packaged 6 per box




Part Number: **72204276**

| | | |
|-------------------|------------------|---------------|
| W 4.3 mm | CT 0.6 mm | TPI 36 |
| UL 35.6 mm | MT 0.4 mm | |



Part Number: **72204228**

| | | |
|-------------------|------------------|---------------|
| W 4.5 mm | CT 0.6 mm | TPI 34 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204287**

| | | |
|------------------|------------------|---------------|
| W 5.8 mm | CT 1.0 mm | TPI 22 |
| UL 8.9 mm | MT 0.6 mm | 90° |

Sagittal Saw Blades - Keyed


▶ Blades are sterile and packaged 6 per box

DISPOSABLES




Part Number: **72204251**

| | | |
|-----------------|------------------|---------------|
| W 5.8 mm | CT 0.6 mm | TPI 22 |
| UL 15 mm | MT 0.4 mm | |




Part Number: **72204229**

| | | |
|-------------------|------------------|---------------|
| W 5.8 mm | CT 0.6 mm | TPI 22 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204248**

| | | |
|-------------------|------------------|---------------|
| W 5.8 mm | CT 0.6 mm | TPI 36 |
| UL 33.5 mm | MT 0.4 mm | |




Part Number: **72204270**

| | | |
|-----------------|------------------|---------------|
| W 9.4 mm | CT 0.6 mm | TPI 36 |
| UL 15 mm | MT 0.4 mm | |



Part Number: **72204235**

| | | |
|-------------------|------------------|---------------|
| W 9.5 mm | CT 0.6 mm | TPI 18 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204271**

| | | |
|-------------------|------------------|---------------|
| W 9.4 mm | CT 0.6 mm | TPI 18 |
| UL 25.4 mm | MT 0.4 mm | |

| | | | | |
|----------|---------------|---------------|--------------------|----------------|
| W | UL | CT | MT | TPI |
| Width | Usable Length | Cut Thickness | Material Thickness | Teeth Per Inch |


Sagittal Saw Blades - Keyed

▶ Blades are sterile and packaged 6 per box




Part Number: **72204226**

| | | |
|-------------------|------------------|---------------|
| W 9.4 mm | CT 0.6 mm | TPI 36 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204284**

| | | |
|-------------------|------------------|---------------|
| W 9.4 mm | CT 0.6 mm | TPI 36 |
| UL 25.4 mm | MT 0.4 mm | (non-sterile) |




Part Number: **72204239**

| | | |
|-----------------|------------------|---------------|
| W 9.4 mm | CT 0.6 mm | TPI 36 |
| UL 41 mm | MT 0.4 mm | |




Part Number: **72204261**

| | | |
|-------------------|------------------|---------------|
| W 9.5 mm | CT 1.0 mm | TPI 36 |
| UL 25.4 mm | MT 0.8 mm | |



Part Number: **72204279**

| | | |
|------------------|------------------|---------------|
| W 10.4 mm | CT 0.5 mm | TPI 36 |
| UL 8.1 mm | MT 0.5 mm | 70° |




Part Number: **72204285**

| | | |
|-----------------|------------------|---------------|
| W 14 mm | CT 1.0 mm | TPI 18 |
| UL 16 mm | MT 0.6 mm | 90° |

DISPOSABLES


Sagittal Saw Blades - Keyed

▶ Blades are sterile and packaged 6 per box



Part Number: **72204244**

W 14 mm CT 0.6 mm TPI 22
 UL 25.4 mm MT 0.4 mm




Part Number: **72204232**

W 14 mm CT 0.6 mm TPI 22
 UL 41 mm MT 0.4 mm



Part Number: **72204245**

W 19.1 mm CT 0.6 mm TPI 22
 UL 25.4 mm MT 0.4 mm



Part Number: **72204241**

W 19.1 mm CT 0.6 mm TPI 22
 UL 41 mm MT 0.4 mm

DISPOSABLES

Reciprocating Saw Blades

▶ Blades are sterile and packaged 6 per box



Part Number: **72204265**

| | | |
|-------------------|------------------|---------------|
| W 6.4 mm | CT 0.6 mm | TPI 20 |
| UL 22.9 mm | MT 0.4 mm | |



Part Number: **72204269**

| | | |
|-------------------|------------------|---------------|
| W 7.4 mm | CT 0.6 mm | TPI 36 |
| UL 23.6 mm | MT 0.4 mm | |



Part Number: **72204246**

| | | |
|-------------------|------------------|---------------|
| W 7.6 mm | CT 0.6 mm | TPI 33 |
| UL 25.4 mm | MT 0.4 mm | |



Part Number: **72204268**

| | | |
|-------------------|------------------|---------------|
| W 8.0 mm | CT 0.7 mm | TPI 32 |
| UL 38.1 mm | MT 0.4 mm | |



Part Number: **72204249**

| | | |
|-------------------|------------------|---------------|
| W 8.0 mm | CT 0.7 mm | TPI 32 |
| UL 25.4 mm | MT 0.4 mm | |




Part Number: **72204256**

| | | |
|-------------------|------------------|---------------|
| W 8.1 mm | CT 0.5 mm | TPI 22 |
| UL 25.4 mm | MT 0.4 mm | |

Reciprocating Saw Blades


▶ Blades are sterile and packaged 6 per box

Part Number: **72204280**



| | | | | | |
|----|-------|----|--------|-----|----|
| W | 9 mm | CT | 0.8 mm | TPI | 24 |
| UL | 40 mm | MT | 0.5 mm | | |

Part Number: **72204250**



| | | | | | |
|----|---------|----|--------|-----|----|
| W | 9.4 mm | CT | 0.6 mm | TPI | 33 |
| UL | 39.1 mm | MT | 0.4 mm | | |



Rasps, Twist Drills, K-Wires, Steinmann Pins & Pin Protectors

Rasps

34 Rasps

Twist Drills

Sterile

35 Synthes

35 Jacobs

Non Sterile

36 Synthes

36 Jacobs

K-Wires

37 K-Wires

Steinmann Pins

38 Steinmann Pins

Pin Protectors

39 K-Wire & Steinmann
Pin Protectors



Rasps - Sterile

▶ Packaged 6 per box

Straight Cut Rasps

| Size | Width | Usable Length | Part Number | |
|-------|--------|---------------|-----------------|--|
| Small | 5.5 mm | 9.5 mm | 72204259 | |
| Large | 6.5 mm | 12.5 mm | 72204267 | |

Cross Cut Rasps

| Size | Width | Usable Length | Part Number | |
|-------|--------|---------------|-----------------|--|
| Small | 5.5 mm | 9.5 mm | 72204233 | |
| Large | 6.5 mm | 12.5 mm | 72204238 | |





Twist Drills - Sterile

▶ Packaged 6 per box

| Sterile Synthes AO-Style Twist Drills | | | |
|---------------------------------------|--------|---------------|-------------|
| Diameter | Length | Usable Length | Part Number |
| 1.5 mm | 85 mm | 60 mm | 72204469 |
| 2.0 mm | 85 mm | 27 mm | 72204454 |
| 2.0 mm | 100 mm | 75 mm | 72204446 |
| 2.0 mm | 125 mm | 100 mm | 72204472 |
| 2.5 mm | 110 mm | 85 mm | 72204443* |
| 2.7 mm | 100 mm | 75 mm | 72204465 |
| 2.7 mm | 125 mm | 100 mm | 72204464 |
| 3.2 mm | 145 mm | 120 mm | 72204444 |
| 3.2 mm | 195 mm | 170 mm | 72204467 |
| 3.5 mm | 110 mm | 85 mm | 72204458 |
| 3.5 mm | 195 mm | 170 mm | 72204460 |



*Titanium nitride coated

| Sterile Jacobs-Style Twist Drills | | | |
|-----------------------------------|--------|-------------|--|
| Diameter | Length | Part Number | |
| 1.6 mm | 127 mm | 72204438 |  |
| 2.0 mm | | 72204431 | |
| 2.0 mm | | 72204445† | |
| 2.4 mm | | 72204434 | |
| 2.8 mm | | 72204435 | |
| 2.8 mm | | 72204456† | |
| 3.2 mm | | 72204433 | |
| 3.2 mm | | 72204450† | |
| 3.6 mm | | 72204470† | |
| 3.6 mm | | 72204455 | |
| 4.0 mm | 127 mm | 72204451 |  |
| 4.4 mm | | 72204462 | |
| 4.8 mm | | 72204457 | |
| 6.4 mm | | 72204468 | |

†European Compression Technique

DISPOSABLES




Twist Drills - Non-Sterile

▶ Packaged 6 per box

Non-Sterile Synthes AO-Style Twist Drills



| Diameter | Length | Usable Length | Part Number |
|----------|--------|---------------|-------------|
| 2.0 mm | 100 mm | 75 mm | 72204442 |
| 2.0 mm | 125 mm | 100 mm | 72204463 |
| 2.5 mm | 180 mm | 155 mm | 72204447* |
| 2.5 mm | 110 mm | 85 mm | 72204439* |
| 2.7 mm | 100 mm | 75 mm | 72204449 |
| 3.2 mm | 145 mm | 120 mm | 72204459 |
| 3.5 mm | 110 mm | 85 mm | 72204448 |
| 4.5 mm | 145 mm | 120 mm | 72204471 |



*Titanium nitride coated

Non-Sterile Jacobs-Style Twist Drills

| Diameter | Length | Part Number |
|----------|--------|-------------|
| 1.6 mm | 127 mm | 72204440 |
| 2.0 mm | | 72204432 |
| 2.4 mm | | 72204436 |
| 2.8 mm | | 72204441 |
| 3.2 mm | | 72204437 |
| 3.6 mm | | 72204452 |
| 4.0 mm | 127 mm | 72204461 |
| 4.4 mm | | 72204466 |
| 4.8 mm | | 72204453 |



K-Wires - Sterile / Non-Sterile

▶ Packaged 6 per box

▶ Organized by trocar points at both ends, one end or threaded

| Sterile K-Wires | | | | | | | |
|-----------------|---------------|----------|----------------|----------------|----------------|----------------|----------------|
| Length | Point Style | Type | Wire Diameter | | | | |
| | | | 0.7 mm / .028" | 0.9 mm / .035" | 1.1 mm / .045" | 1.4 mm / .054" | 1.6 mm / .062" |
| 4" | Single Trocar | Smooth | | 72204558 | 72204549 | | 72204547 |
| | Double Trocar | Smooth | 72204534 | 72204524 | 72204510 | 72204535 | 72204515 |
| 6" | Double Trocar | Smooth | 72204537 | 72204522 | 72204507 | 72204530 | 72204509 |
| | | Threaded | | | 72204557 | | 72204548 |
| 9" | Single Trocar | Smooth | | | 72204540 | 72204541 | 72204529 |
| | Double Trocar | Smooth | | 72204538 | 72204518 | 72204552 | 72204514 |
| | | Threaded | | | 72204556 | | 72204543 |

| Non-Sterile K-Wires | | | | | | | |
|---------------------|---------------|----------|----------------|----------------|----------------|----------------|----------------|
| Length | Point Style | Type | Wire Diameter | | | | |
| | | | 0.7 mm / .028" | 0.9 mm / .035" | 1.1 mm / .045" | 1.4 mm / .054" | 1.6 mm / .062" |
| 4" | Single Trocar | Smooth | 72204555 | 72204542 | 72204523 | 72204533 | 72204526 |
| | Double Trocar | Smooth | 72204527 | 72204520 | 72204508 | 72204525 | 72204512 |
| 5.5" | Single Trocar | Smooth | | 72204554 | 72204545 | | 72204550 |
| 6" | Double Trocar | Smooth | 72204546 | 72204532 | 72204517 | 72204536 | 72204521 |
| 9" | Double Trocar | Smooth | | 72204528 | 72204513 | 72204551 | 72204511 |
| | | Threaded | | | 72204553 | | 72204544 |
| | Single Trocar | Smooth | | 72204539 | 72204531 | | 72204519 |

DISPOSABLES



Steinmann Pins - Sterile / Non-Sterile

▶ Packaged 6 per box

▶ Organized by trocar points at both ends, one end or threaded


Sterile Steinmann Pins

| Length | Point Style | Type | Wire Diameter | | | | | | |
|--------|----------------|--------|----------------|----------------|----------------|---------------|----------------|----------------|----------------|
| | | | 2.0 mm / 5/64" | 2.4 mm / 3/32" | 2.8 mm / 7/64" | 3.2 mm / 1/8" | 3.6 mm / 9/64" | 4.0 mm / 5/32" | 4.8 mm / 3/16" |
| 9" | Single Trocar | Smooth | 72204570 | 72204561 | 72204568 | 72204577 | 72204585 | 72204576 | 72204588 |
| | Double Trocar | Smooth | 72204579 | 72204587 | | | | | |
| | Single Diamond | Smooth | 72204580 | | | | | | |
| | Double Diamond | Smooth | 72204582 | | | | | | |

Non-Sterile Steinmann Pins

| Length | Point Style | Type | Wire Diameter | | | | | | |
|--------|----------------|----------|----------------|----------------|----------------|---------------|----------------|----------------|----------------|
| | | | 2.0 mm / 5/64" | 2.4 mm / 3/32" | 2.8 mm / 7/64" | 3.2 mm / 1/8" | 3.6 mm / 9/64" | 4.0 mm / 5/32" | 4.8 mm / 3/16" |
| 9" | Single Trocar | Smooth | 72204560 | 72204559 | 72204569 | 72204564 | 72204575 | 72204578 | 72204574 |
| | Double Trocar | Smooth | 72204562 | 72204565 | 72204571 | 72204583 | | | |
| | | Threaded | 72204584 | | | | | | |
| | Single Diamond | Smooth | 72204566 | 72204563 | 72204573 | 72204567 | 72204581 | | |
| | Double Diamond | Smooth | 72204572 | 72204586 | | | | | |

DISPOSABLES



Pin Protectors - K-Wire and Steinmann Pin

▶ 4 per package, 6 packages per box


K-Wire Pin Protectors

| | 0.9 mm / .035" | 1.1 mm / .045" | 1.6 mm / .062" |
|--------------|-----------------|-----------------|-----------------|
| | 72204499 | 72204497 | 72204498 |
| Color | Blue | White | Green |

Steinmann Pin Protectors

| | Wire Diameter | | | | | | |
|--------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| | 2.0 mm / 5/64" | 2.4 mm / 3/32" | 2.8 mm / 7/64" | 3.2 mm / 1/8" | 3.6 mm / 9/64" | 4.0 mm / 5/32" | 4.8 mm / 3/16" |
| | 72204501 | 72204502 | 72204506 | 72204503 | 72204505 | 72204504 | 72204500 |
| Color | Pink | Grey | Tan | Lavender | Brown | Navy Blue | Orange |

DISPOSABLES



Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service: +1 800 343 5717
International Customer Service: +1 978 749 1140

®Trademark of Smith & Nephew. Certain Marks Registered US Patent and Trademark Office. All trademarks acknowledged. ©2014 Smith & Nephew. All rights reserved. Printed in USA. 02185 V1 07/14