

EXCEL G-LOGIC Specification Sheet

Equipped with Attendant brakes

Interchangeable upholstery colours

Half-folding backrest

Height adjustable armrests

Seatbelt as standard

Swing-away and quick release legrests

Quick release wheels for transportation


4 frame colours available

Anti-tip wheels as standard

Puncture proof tyres



G-Logic 



Easy to remove
Attendant brakes



Anti-tip wheels



Transit version
available

Available Frame Colours



Upholstery Colours



Full Specifications

Seat Width	40, 45, 50 cm
Seat Depth	43 cm
Seat Height	50 cm
Backrest Height	44 cm
Total Length	95 - 104 cm*
Total Width	62, 67, 72 cm
Total Height	95 cm

Total Weight	12 - 15 kg**
Maximum User Weight	113 kg - 135 kg**
Transit Available	Yes
Self-propel Available	Yes
Crash Tested	Yes

KEY: * = Adjustable or optional ** = Varies depending on seat width: 113 kg (17 Stone) for 40 & 45 cm seat width, 135 kg (21 Stone) for 50 cm seat width



www.vanosmedical.co.uk