

Specifications	
Chair	
Weight Capacity	18 st 8 lb (120 kg)
Overall Length	1100 mm (43.3")
Overall Width	620 mm (24.4")
Drive Tyres	Front: 195 mm (7.6") Rear: 250 mm (9.8")
Ground Clearance	40 mm
Suspension	Limited
Colour(s)	Silver & Black
Power	
Turning Radius	650 mm (25.59")
Max Speed	5 mph (8 km/h)
Drivetrain	2 x 24v / 200w Brushed Motors
Dual Breaking System	Electromechanical
Battery	
Battery Requirements	1 x 20 Ah Lithium Battery
Battery Weight	4.6 kg
Battery Charger	5 - Amp, Off-Board
Total Weight Without Batteries	27.8 kg

Features



Compact
Lightweight Design



Remote Controlled Automatic
Folding System



Back Seat Pocket For
Extra Storage While on-the-go

TECHNICAL SPECIFICATIONS

MODEL:SW6000E

Item		Standard Requirements
Max.Speed		8.5KM/H
Braking Performanc	Horizontal road brake	$\leq 1.5M$
	Maximum safe slope braking	$\leq 2M(6^\circ)$
Slope performance		12°
Static stability	Forward leaning	$\geq 9^\circ$
	Backward leaning	$\geq 9^\circ$
	Side leaning	$\geq 9^\circ$
Dynamic stability	Backward dynamic stability of the ramp	$\geq 6^\circ$
	Forward dynamic stability of the ramp	$\geq 6^\circ$
	Side dynamic stability of the ramp	$\geq 6^\circ$
Ground Clearance		$\geq 40MM$
Wedge width		100MM
Gradient		9°
Turning radius		650MM
Driving Range		$\geq 20KM$
Length*Width*Height		1100*620*960MM
Length*Width*Height(Folded size)		630*355*770MM
Front wheel		195*45MM
Rear wheel		250*60MM
Motor		Brush 24V/200W
Battery		Lithium battery

REVIEWER: Kimi Zhang

DATE: 20th,Jan,2019

COMPANY STAMP:

地址: 中国·苏州市·工业园区新发路 68 号

Add:Building No.68 XinFa Road,Suzhou Industrial Zone,Suzhou.

电话: 0512-82283862

Tel:0512-82283862

传真: 0512-82283866

Fax:0512-82283866

邮编: 215123

Postal code:215123