## KING COBRA



Following the success of our Cobra scooter we've developed the new King Cobra. This offers an evolved dynamic design, generous spaciousness, a 31 mile range and a maximum weight limit of 32 stone.



The King Cobra is a revolutionary scooter in terms of design, aesthetics and styling. It has been aerodynamically designed to give an outstandingly stable drive and maximise battery consumption. It's a remarkable scooter that makes a huge statement and sets a new benchmark with its striking dynamic design.



# HE SET





All round suspension

# Superb Comfort, Style & Performance...

- » Top speed of 8mph\*.
- » Maximum range of 32 miles on full battery charge\*.
- » Stylish & distinctive design.
- » Comfortable rotating seat with sliding and reclining facility.
- » Height adjustable headrest and angle adjustable armrests to suit user requirements.
- » Electromagnetic brake system and emergency cable brakes.
- » All round suspension ensures a smooth ride and optimal handling.
- » Generous tiller adjustment to suit user requirements and assist easy transportation.
- » User friendly controls with a modern digital LCD display.
- » Automatic speed reduction when steering for added safety.
- » Twin USB power socket.
- » Fully functioning indicators and twin adjustable rear view mirrors for extra visibility.

- » High level indicators for added safety.
- » Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark.
- » Stylish 'Blade' alloy wheels with black pneumatic tyres for a smoother ride.
- » Delta bars as standard.
- » Bottle/cup holder included.
- » Rear anti tip wheels as standard.
- » Anti roll back system safety device prevents scooter from rolling backwards on hills.
- » Freewheel facility allows movement of scooter without switching the motor on.
- » Optional rear pannier box.
- » Recommended MK Gel batteries: 2 x 12V 96Ah.
- » Weight capacity 200kg (32st).
- Subject to conditions of use.

### CODE: KCOBRAGG DESCRIPTION: King Cobra



Length: 170cm (67")



Height: 140cm (54")



Width: 73cm (29")



Max Capacity: 200kg (32st)



Ground Clearance: 11.5cm (4.5")



Max Gradient:



Stylish 'Blade' Alloy Wheels



LED Lighting

### Leaflet Code: LL208

If you would like to view this product please contact your local stockist, details below.