

## Woodpecker Endo Radar Plus Programmes – SGP Mode

Order	Manufacturer	System
1	FKG	Endo Shaper
2	FKG	Finisher
3	FKG	Gates
4	FKG	BioRace
5	FKG	iRace
6	Micro Mega	Endo Flare
7	Micro Mega	G-Files
8	Micro Mega	One Shape Apical
9	Micro Mega	One Shape
10	Micro Mega	Revo-S
11	SybronEndo	TF
12	SybronEndo	K3XF
13	SybronEndo	K3
22	Reciprocating	Safety Glide Path
23	Program	Open

Includes - 8x Open Programmes      Fully Programmable Speed/Torque/File/Size

Plus

- (A) 1x SGP Open Programme (Reciprocating)
- (B) Reciprocating  
Angles x3  
Adjustable from 20°- 400° in 10 ° intervals.