

Brand	File system	File No.	Speed	Torque	Forward angle	Reverse angle	
WOODPECKER	W3_Pro	17/12	300	3.0			
		18/05	350	2.0			
		25/06					
		12/vt					
		16/vt	350	1.5			
		35/vt					
	W3_ONE	20/07				30	150
		25/08	400	4.0			
		30/09					
	W3_Single	S1					
		S2	500	2.0			
		S3					
	W2_Plus	N0					
		N1	350	3.0			
		W2					
		W3	350	2.0			
		W1					
	W2_Pro	N2	350	1.5			
		17/12	300	3.0			
		20/04					
	W2_One	25/06					
		35/04	350	2.0			
		16/02					
		19/02					
		20/06					
	W2	25/08	400	4.0			
		40/08					
		17/10					
W2	20/04						
	25/06						
	35/04	450	2.5				
	17/10						
	20/04						
	25/06						
Endo_Plus	35/04						
	17/12	300	3.0				
	20/06	350	2.0				
	25/06						
	20/04	350	1.5				
Dentsply	18/05	350	2.0				
	12/vt	350	1.5				
	WAVE ONE	S				30	150
		P	400	4.0			
		L					
Wave_One_Gold	20/07				30	150	
	25/07	400	4.0				
	35/06						
	45/05						
PROTAPER	S1	250	3.0				
	SX						
	D3	250	1.5				
	D1						
	D2						
	F3						
	F4	250	2.0				
	F5						
	F2						
	F1	250	1.5				
	S2	250	1.0				
PROTAPER_NEXT	P1						
	P2						
	X4						
	X5	300	2.0				
	X1						
	X2						
PROTAPER GOLD	X3						
	SX						
	S1						
	S2						
	F1	300	3.0				
	F2						
	F3						
VDW	F4						
	F5						
	F25						
	R25						
	R40	400	4.0	50°	170°		
RECIPROC	R50						
	25/08						
	40/06	400	4.0	50°	170°		
RECIPROC BLUE	50/05						
	10/04						
Mtwo	30/05	250	1.5				
	30/06						
	35/06						
	60/04						
	25/06	250	3.0				
	25/07						
	50/04						
	45/04						
	20/06	250	2.5				
	40/04						
	40/06						
VDW	35/04	250	2.0				
	15/05						
	70/02						
	40/06	250	3.0				
	50/02						
	60/02						
VDW	35/06	250	2.5				
	45/02						

Brand	File system	File No.	Speed	Torque	Forward angle	Reverse angle	
SybronEndo	K3	15/02					
		20/02	350	1.0			
		25/02					
		30/02					
		35/02	350	1.5			
		15/04					
		40/02	350	2.0			
		20/04					
		25/04					
		30/04					
		35/04					
		15/06					
		20/06	350	2.5			
		25/06					
		30/06					
	35/06						
	45/02						
	40/04						
	45/04						
	50/04						
	55/04						
	60/04						
	40/06	350	3.0				
	45/06						
	50/06						
	55/06						
	60/06						
	25/08						
25/10							
25/12							
SybronEndo	K3XF	15/04					
		20/04					
		25/04					
		30/04	350	2.0			
		35/04					
		15/06					
		20/06					
		25/06					
		30/06					
		35/06					
		40/04					
		45/04					
		50/04					
		55/04					
		60/04					
	40/06	350	3.0				
	45/06						
	50/06						
	55/06						
	60/06						
	20/08						
	25/08						
	30/08						
	40/08						
	20/10						
	25/10						
	30/10						
	40/10						
20/12							
25/12							
30/12							
40/12							
MICRO-MEGA	Revo-S	25/04	500	2.5			
		40/04					
		25/06					
		30/06					
		25/08	500	3.0			
MICRO-MEGA	One Shape	25/10					
		25/12					
		SC1					
		SC2					
		SU	400	1.5			
MICRO-MEGA	One Shape Apical	AS30					
		AS35					
		AS40					
		25/06	400	2.5			
		Apical2					
MICRO-MEGA	G-Files	Apical1	400	1.5			
		G1					
		G2	400	2.0			
		25/12	400	2.5			
		EndoFlare					
MICRO-MEGA	iRace	R1					
		R2					
		R3	600	1.5			
		R1a					
		R1b					
MICRO-MEGA	iRace	BR0					
		BR1					

	FlexMaster	30/06	250	2.0		
		40/04				
		35/04				
		40/02	250	1.5		
		25/06				
		30/04				
		35/02	250	1.0		
		20/06				
		25/04				
		30/02				
		20/02				
		20/04				
		25/02				
		15/02				
		15/04				
15/06						

FKG	BioRace	BR2	600	1.0		
		BR3				
		BR4				
		BR5				
	Gates	1	650	1.0		
		2				
		3	650	3.0		
		4				
		5				
6						
XP-endo Finisher	25/00	800	1.0			
XP-endo Shaper	30/04	800	1.0			
User defined continuously rotate system	PROGRAM	1	250	1.0		
		2	250	1.5		
		3	250	2.0		
		4	250	2.5		
		5	250	3.0		
		6	250	3.5		
		7	250	4.0		
		8	600	0.8		
User defined reciprocation system	RECIPROCATING	1	400	4.0	30	150
		2	400	4.0	50	170
		3	400	4.0	60	60