

TABURETT Assistant

Specially developed to be used by assistants, the armrest provides flexibility and reach thanks to improved support to back, arms and abdominal areas. A high gas spring and foot ring offers easy adjustability, and the deep seat can be tilted to suit your needs.

Description

Seat height	56-91 cm
Seat dimensions	37x44 cm
Gas spring options	3
Alternative castors	Yes
Base diameter	58 cm
Adjustable seat	Yes, 15°
Upholstery	Scandinavian leather, vinyl or fabric
Weight	16 kg
Max. load	120 kg
Warranty	5 years

Accessories

Foot control	Foot operated height adjustment
--------------	---------------------------------

Seat Heights

User height	Gas spring	Seat height
< 180 cm	Standard	56-65 cm
180-190 cm	High	63-77 cm
> 190 cm	Extra High	71-91 cm



Key Features

Mechanism

Two lever spring adjusted mechanism allows for adjustable height and up to 15° lateral movement.

Seat

Made with 100% recyclable, injection moulded plastic (PA6).

Padding

Upholstered with flexible polyurethane foam, recyclable when using a certain process.

Armrest padding

Upholstered with injection moulded foam, recyclable when using a certain process.

Upholstery

Our standard vinyl contains 91.5% PVC and 8.5% Polyamide and is fire resistant according to class 1/M-M2. The leather is high quality Scandinavian aniline or semi-aniline leather. The manufacture is certified to environmental standard SS-EN ISO 14001-1996.

Gas spring

Columns are available in three sizes. The spring is made of polished steel in Germany using nitrogen, a sustainable gas existing in 78% of normal air. Gas springs are maintenance-free and recyclable.

Armrest

Support arm is made from wood, mechanism from metal.

Base

Steel.

Castors

Cores are made from polyamide 6 and tread is made of recyclable polyurethane.