



TURBO-JET COMPACT

Control the power

Cattani's Turbo-Jet Compact is an extremely robust suction system, for single surgery dentists requiring a fixed speed, reliable and cost-effective solution. Although small and light, this machine is easy to use and delivers a powerful performance.

This unit is ideal for general conservation dentistry, providing very high airflow and a vacuum power capable of lifting debris with a larger bore tip.

Supplied with built-in amalgam separation and pre-wired for installation, this suction system easily fits into a standard surgery cabinet. Cattani Turbo-Jet Compact comes with a one year guarantee and Cattani's vacuum pump protection.

Designed for single surgery use

Call us: 01527 877997

e: info@cattaniesam.co.uk

w: www.cattaniesam.co.uk/smart





TURBO-JET COMPACT



Fixed speed suction and very high airflow

- 4 amp current and rated 0.42 kW
- Consistent suction at 650 l/min and 160 mbar vacuum
- 61 dB(A) noise level, reducing to 52 dB(A) with sound proof box
- Extremely reliable, durable and proven product



Built-in amalgam separator

- Up to 95.6% retention rate at 5.5 l/min maximum
- Complies to ISO 11143 – Dental Equipment Amalgam Separators
- Amalgam container is easy to access and change
- Optional bacteria/mercury filter for further cross-infection control



Simple and fast installation and maintenance

- Supplied pre-wired for installation
- Fits in 300mm cabinet
- Dimensions: H 495mm / D 390mm / W 250mm / weight 33 kgs
- Easy to access, reliable components

Supplied with a two year warranty

- Designed for longevity
- Motor protected by thermal device
- Cattani vacuum pump protection
- Rapid and reliable technical support and after sales service

For further information, please contact your local stockist: