

# Murray's

Medical Seating

Winners of the Clinical Seating OJEU contract  
for Barts & The London NHS  
Criteria: Ergonomics, Infection Control and  
Build Quality

Operating Theatres  
Anaesthetics  
Eye Clinics  
Recovery  
ICU – HDU  
Pain Management  
Ultra sound  
Laboratories  
Maternity Units  
Delivery Suites  
Maxillofacial  
Dental  
All Clinical Areas



**Colin Sole** General Manager  
• T: 01243 811881 • E: colin@murrayequipment.co.uk



British Made



British Made


# Murray's Medical Seating

► Murray Equipment Co. Ltd. Was established in 1953 and is a respected manufacturer of medical seating for clinical departments, which supplies directly to the NHS and dental clinics.

► Murrays provide exceptional quality, and great value. We are so confident with the quality of our products, we offer a 5 year warranty\*. This in turn results in large savings, as you will not have to replace them as regularly as other stools. Furthermore, our ergonomically designed seating will reduce back pain, fatigue, and long term RSI issues, which will save vast amounts of money over the years by reducing absences from highly paid and trained clinicians. This was highlighted by a three year clinical study at Nottingham University Hospital; this study was only completed with one clinical therapy area, yet applies to all clinical therapy areas, whereby we manufacture specific chairs for those areas.

► All products comply with fire safety regulations, CE Medical Devices 93/42/EEC and ISO 9001. All of our seating is free from latex. The vinyl we use in impregnated with anti bacterial agents during the manufacturing process and an anti fungal backing is used on the material. Our wipe clean bases are strengthened by a T bar which makes it very easy to clean, compared to a typical plastic honeycomb base which is perfect place for bacteria to breed. Murray bases are chrome, so that any dirt, dust or bodily fluids are easily identifiable for cleaning purposes.

► Murray's seating is handmade at our factory in Arundel, therefore all our chairs are made to order. A range of optional extras and a wide choice of colours means that you have the flexibility to purchase the seating you need to meet your exact specification. Colours are often used to identify which department or surgery the chair belongs to.



\*Warranty does not apply to misuse or upholstery. Also Luna only has a 2 year warranty.

# Murray's Medical Seating



British Made





# Helping you choose the right seating

Throughout the brochure you may come across terms that do not mean anything to you, therefore please see our jargon buster below, which is intended to assist you to choose the correct seating. The use for each model that is described throughout the brochure is only intended as a guide.

## ► Height Range

This is the adjustable height range of the cylinder from the lowest to the highest point, measured from the floor to the top of the seat pad.

Our seating is normally supplied with the option of a standard or high cylinder; however, if required, we are able to supply a lower or higher option. Should this be required, please contact us for details. An adjustable foot ring/rest is supplied as standard when ordered with a high cylinder - whilst not recommended this may be removed if required.



High cylinder highest point



Standard cylinder highest point



High cylinder lowest point



Standard cylinder lowest point

Our chairs are fully adjustable so please do not dismiss the comfort of a chair until it's been fully adjusted to suit the individual user.

## ► Weight capacity - Gas cylinder

Cylinders have been tested to hold a weight of 22 stones but this may be increased if required.

Detail are available on request.

## ► Tilt facility

The seat pad tilts back and forth to enable pelvic movement, which helps promote a good neutral spinal position.

On all models, where applicable, the back rest and the seat pad can be independently adjusted and used in a free floating position or locked.





# Helping you choose the right seating

## Fully adjustable Back Rest

### ► Rake Adjustment

Spring loaded adjustable movement. This can be in a free floating position, which means, as the operative leans back and forth, the back rest will move with them providing continued support; alternatively it can be locked into position.



### ► Height Adjustable Backrest

This is raised manually up and down, to achieve the correct position for the lumbar support, and then locked into place



### ► Double Curvature

This means a contoured moulded seat pad.

### ► Restraint Castor

These act on vinyl flooring, similar to a normal free wheeling castor would act on a carpeted floor. This helps prevent the chair/stool skidding away from the operative, whilst giving freedom of movement.



## ► Foot Height Control (FHC)

This is an optional extra (see pages 20-21) that enables the height range of the seat to be adjusted with the foot, instead of a hand lever, ideal use for when in scrubs. This can be fitted to all models.



## ► Material

All products comply with fire safety regulations and ISO 9001. All our seating is free from latex. The upholstery is breathable, wipe clean and impregnated with antibacterial agents at the point of manufacture. Murrays offer a choice of 85 colours, all included in the price (see page 22). It is common practice that the colour is used as a device to identify which department or theatre the chair/stool belongs to.

## ► The Base

of the chair is polished aluminium.

Many other bases are strengthened by a honeycomb effect on the underside creating an ideal place for bacteria to breed.



Polished aluminium base

Our base is strengthened by a T structure which gives a smooth clean surface to wipe clean.

(The hydraulic range and 19" base has a very basic honeycomb effect to strengthen it).



Typical honeycomb base

All our chairs/stools are fitted as standard with **Restraint Castors or Static feet**. However, we have a full range of locking castors; please refer to the optional extras on pages 20-21 for a detailed description of their functions.

Trials available for two week test periods. A small charge will be applied. Please contact Colin at Murrays for more information.





HCAP-GT

- ▶ Furthermore the back rest has an incorporated ratchet back system, meaning it can be adjusted using forearms, so if wearing gloves or scrubs you do not have to touch the back.
- ▶ The chair has a plastic casing covering the back and the underside, allowing for protection against bacteria and germs and an easy wipe surface.
- ▶ The chair features removable seat and back pads for further ease of cleaning and a cost effective way of replacing damaged/dirty pads.
- ▶ Covers many clinical areas.

The Capricorn is a multi functional seat, viable for use throughout a hospital, dental clinic or laboratory.

- ▶ The chair features a double curvature seat pad, with injected moulded foam, providing a comfortable working position for long periods of time.
- ▶ It has a fully adjustable backrest with tilt function; there is also a tilt function in the seat pad allowing for a versatile range of movement to find each specialists specific sedentary position.
- ▶ The Capricorn is sold at 2 different working heights which can be discussed over the phone, to best suit your needs.



SCAP-GT

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.



HCAP-GT

Code	Model Description
HCAP-GT	Capricorn Chair: Adjustable tilt and backrest. Adjustable height cylinder and foot ring.
Height range	520 - 720mm Approximate height ranges are measured from floor to top of seat pad.



SCAP-GT

Code	Model Description
SCAP-GT	Capricorn Chair: Adjustable tilt and backrest. Adjustable height cylinder.
Height range	520 - 720mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



SCAP-GT - Colour variations



British Made

# Gemini Range



SGEM-GT



Adjustable arms

- see optional extras on page 21

360° swing around backrest



► The SGEM-GT. This is a traditionally shaped seat and is a little wider than the other ranges. It has full independent adjustment on the seat tilt and back rake. This chair can be fitted with most of the optional extras on pages 20-21. Ideal for Anaesthetists, Surgeons, Dentists and Ophthalmologists.



ratchet system

► The HGEM-SAB is a specialist stool. Whilst the back rest will act as a lumbar support, it has been designed for an operative that requires a moveable arm, without having to leave the stool.

The arm can rotate 360 degrees, which provides an arm with one swift movement. Also the arm has a ratchet system allowing it to move closer or further away from the operative. It will also wrap a round in front of the operative providing a torso butt, without having to straddle the seat.

This is very popular with Dental Nurses, Head & Neck Surgeons and Vascular Sonographers.



HGEM-SAB


Round stool with 360° swing around backrest.

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.






# Gemini Range




HGEM-SAB

Code	Model Description
HGEM-SAB	Round stool with 360° swing around backrest. Adjustable high cylinder and foot ring. No tilt.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.




SGEM-SAB

Code	Model Description
SGEM-SAB	Round stool with 360° swing around backrest. No tilt. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.



HGEM-GT

Code	Model Description
HGEM-GT	Gemini chair. Adjustable tilt and backrest. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.



SGEM-GT

Code	Model Description
SGEM-GT	Gemini chair. Adjustable tilt and backrest. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



# Advance Range



HADV-GN

► The HADV-GN is a round stool with a fully adjustable back rest, without seat pad tilt facility, and is available with standard or high cylinder.

Ideal for Laboratories, Maternity, Eye Clinics, Theatre Recovery and as a general clinical stool.

This chair works really well for Anaesthetists when it is fitted with the table attachment - see optional extras on page 21.

Two-tone colour option available at no extra cost on all chairs.

This is a multi functional range of seating, providing a solution for most needs within theatres and other clinical areas.

► The SADV-GT has a double curvature seat pad, with injected moulded foam, providing extreme comfort for the clinician who has to sit for a long period of time.

It has a fully adjustable backrest with tilt, as does the seat pad. The back rest and the seat pad both tilt separately and can be fully adjusted to provide maximum comfort and support. This range can be supplied with a standard or high cylinder. Available with arms as optional extras - see pages 20-21.

Great for Surgeons, Anaesthetists, Dentists and Sonographers.




tilt facility

SADV-GT

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.




# Advance Range




SADV-GT-H

Code	Model Description
SADV-GT-H	DC seat. Adjustable tilt and backrest. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.




SADV -GT

Code	Model Description
SADV-GT	DC seat. Adjustable tilt and backrest. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.



HADV-GN

Code	Model Description
HADV-GN	Round seat with backrest. No tilt. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.



HADV-GN-S

Code	Model Description
HADV-GN-S	Round seat with backrest. No tilt. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.

# Advance Range



British Made



HADV-NB



SADV-NB

Code	Model Description
HADV-NB	Round seat without backrest. No tilt. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.

Code	Model Description
SADV-NB	Round seat without backrest. No tilt. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.





# Saddle Seats


► The Saddle stool helps to maintain an upright posture when seated with the pelvis held securely in a neutral position. The Saddle stool is available with a standard height or high cylinder. Typically used in Radiography, Dentistry, Theatres and in Laboratory departments.



MUR-SAD




MUR-SAD-S



MUR-SAD

Code	Model Description
MUR-SAD	Saddle seat without backrest. Adjustable tilt facility. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.



MUR-SAD-S

Code	Model Description
MUR-SAD-S	Saddle seat without backrest. Adjustable tilt facility. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



# Posture Range

The Posture range has been designed for the operative who needs to work close to the patient. The front edge of the seat has been cut away to encourage the legs to drop sideways rather reaching forward. Whilst this range is likened to the saddle seat, it is not as harsh, therefore providing a broader appeal.



POS-DC



POS-DC  
with Dixon arms & FHC  
- see optional extras on  
page 21



POS-DC-S


► The POS-DC has a fully adjustable back rest and seat tilt facility, enabling a good neutral sitting position.

This is a No. 1 best seller for Sonographers, however is widely used in all types of Surgery, ICU-HDU, Delivery Suits, Eye Clinics and Laboratories when fitted with the Dixon arms and FHC as shown above. It is ideal for Ophthalmic and Micro Surgeons - see optional extras on page 21.

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.



# Posture Range



POS-DC

Code	Model Description
POS-DC	DC seat. Adjustable tilt and backrest. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.



POS-DC-S

Code	Model Description
POS-DC-S	DC seat. Adjustable tilt and backrest. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.



POS-NB

Code	Model Description
POS-NB	DC seat without backrest. Static tilt. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



British Made

# Five Star Range



SLS-DC

► The Five Star range has a double curvature seat pad, which is longer than that of the rest of our range. This gives greater support under the thigh. It has exactly the same chassis as the Advanced GT which provides a fully adjustable back rest and seat tilt facility.

This range will accommodate the Dixon arms - see optional extras on page 21, and is popular with Ophthalmic, Head & Neck and Micro surgeons.



HFS-DC

► HFS-DC This chair has an adjustable high cylinder and foot ring. It is ideally suited to accommodate any of our arms, and is perfect for the fully adjustable, removable Stainless Steel arms (Dix-arm) - see optional extras on page 21, which are used by many ophthalmic, Head & Neck surgeons. When used without arms, it is ideal for Anaesthetists, Sonographers, and Lab staff.

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.





British Made

# Five Star Range



HFS-DC

Code	Model Description
HFS-DC	DC seat. Adjustable tilt and backrest. Adjustable high cylinder and foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.



SLS-DC

Code	Model Description
SLS-DC	DC seat. Adjustable tilt and backrest. Adjustable standard height cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



# Hydraulic

## ▶ CHROMA-HYD PATIENT CHAIR

This model is designed for the patient. It has a large seat pad and back rest and is fitted with secure arms. The arms are shortened and shaped downwards to enable the chair to be pushed under a table, such as a slit lamp, without the patient's hands or arms being trapped under the table. The chair can be lowered to enable the patient to be seated easily and safely, then once the patient is safely seated, the clinician can raise the patient to the required position without the risk of personal injury. The patient can then be lowered to enable a safe exit.

This chair has proved very popular with Eye Clinics for use with the slit lamp. However it has also been used in other clinical departments where the patient needs to be raised and lowered by a clinician.

- ▶ Raise and lower seat height whilst patient seated
- ▶ Ability to lock the swivel action
- ▶ Large seat and backrest for patient use
- ▶ Fixed foot ring/rest
- ▶ Fitted with static feet or castors
- ▶ Height range 490mm to 650mm
- ▶ Weight capacity 32 stones/203kg
- ▶ Delivered fully assembled



CHROMA-HYD



CHROMA-HYD PATIENT CHAIR

Code	Model Description
CHROMA-HYD	Adjustable seat height. Lockable swivel action. Adjustable hydraulic cylinder.
Height range	490 - 650mm Approximate height ranges are measured from floor to top of seat pad.

**HOW TO ORDER** Choose the code for the model required. If any optional extras are needed, add the code for the optional extra/s from pages 20-21. Hyphenate each extra code. For example, if you want to order a chair from our Posture range with load locking castors and FHC then the order would read: POS-DC-LL Castors-FHC. Don't forget to add the colour of your choice - see page 22.

# Luna Range



British Made



S-LUN-GT



H-LUN-GT



H-LUN-GT



S-LUN-GT

Code	Model Description
H-LUN-GT	Adjustable Seat Tilt. Standard round flat seat. Adjustable high cylinder with foot ring.
Height range	High cylinder 540 - 725mm Approximate height ranges are measured from floor to top of seat pad.

Code	Model Description
S-LUN-GT	Adjustable Seat Tilt. Waterfall front. Adjustable standard cylinder.
Height range	Standard cylinder 480 - 625mm Approximate height ranges are measured from floor to top of seat pad.

Discounts available on volume orders. Delivery charge £10.00 per item. UK mainland excluding Scottish Highlands and offshore islands (out of area price on application). Murray stools are covered by our 5 Year Guarantee, excluding upholstery.



# Optional Extras

Our standard chairs are supplied with the restraint castor - see page 5. Glides can be supplied at no extra cost - see below.



Code	Description
LL Castors	<b>Load Locking Castors</b> provide ease of movement from one place to another. Castors lock automatically when weight applied.
un-LL Castors	<b>Un-Load Locking Castors</b> provide stability at the point of sitting in the chair. Castors unlock automatically when weight applied providing movement while working.



Code	Description
65mm-manloc	<b>65mm Manual Locking Castors</b> provide the user with the choice to lock the chair in position or have the castor freewheeling.



Code	Description
Glides	<b>Glides (Static Feet)</b> plastic feet, to provide a static stool at all times, typically used in Laboratories. <b>Supplied in lieu of standard castors.</b>



Code	Description
FHC	<b>Foot Height Control (FHC)</b> This is an optional extra that enables the height range of the seat to be adjusted with the foot instead of a hand lever. Ideal for use when in scrubs. Please note: a high cylinder is used when a FHC is fitted.





British Made

# Optional Extras



Code	Description
Dix-arm	<b>Dixon Arm</b> Stainless steel fully adjustable arm. Precision engineered, laser cut for accuracy, enabling a 360° rotating action at 2 points. This enables the arm to extend forwards or wrap around in front of the operative to provide an elbow rest. The arm is height adjustable, with a good range, and is removable. The arm provides total flexibility. Available individually or as a pair. Available for Five Star, Posture and Gemini.



Code	Description
Adjarm	<b>Adjustable Arm</b> A basic arm that is height adjustable via a ratchet system, as an upgrade and can also be provided with a swivel action. Upholstered to match the chair. Available for Five Star, Posture and Gemini.
Adjarm-SA	Height adjustable arm no swivel
Removable-SA	Height adjustable arm with swivel
	Removable adjustable arm with or without swivel



Code	Description
APC 19	<b>19" Base</b> small base for when space is limited or the clinician needs to get close to the patient.



Code	Description
Ft-ring	<b>Adjustable Foot Ring</b> Supplied as standard with high cylinder. Can be supplied as an optional extra on a standard height cylinder.



# Upholstery Shades

							
Adobe White US 394	Ivory US 500	Coral Sand US 502	Froth US 501	Surf US 513	Fjord Blue US 517	Bristol Blue US 515	Pink US 503
							
Alabaster US 396	Parchment US 385	Shell US 402	Sand US 384	Dusty Jade US 417	Turquoise US 419	Dutch Blue US 518	Mauve US 415
							
Taupe US 410	Bone US 404	Peach US 406	Vicuna US 401	Northwoods US 528	Grotto US 426	Wedgewood US 414	Tea Rose US 416
							
Putty US 403	Bisque US 409	Rose US 412	Doe US 519	Lilac US 508	Sage US 526	Space Blue US 353	Rouge Red US 505
							
Gray Moss US 525	Dewberry US 424	Sunset US 521	Camel US 501	Crocus US 433	Artichoke US 527	Royal US 354	Raspberry US 362
							
Dove US 413	Woodrose US 423	Paprika US 370	Oak US 421	Deep Violet US 511	Olive Green US 529	Regimental Blue US 345	Dark Cherry US 507
							
Greystone US 425	Chestnut US 375	British Tan US 365	Terra Cotta US 376	Grape US 509	Alpine US 530	Imperial Blue US 432	American Beauty US 360
							
Gunmetal US 392	Bronze US 523	Gingersnap US 369	Cinnamon US 377	Sangria US 506	Forest US 344	Sky Blue US 516	Tomato US 358
							
Blue Ridge US 427	Mocha US 381	Cinnabar US 520	Goldenrod US 324	Plum US 361	Yew Tree US 350	Capri Blue US 514	Rosalee US 504
							
Graphite US 429	Rustic Brown US 522	Claret US 357	Golden Corn US 524	Dark Purple US 512	Emerald US 343	China Green US 434	Mandarin Orange US 372
							
Black US 393	Espresso US 383	Burgundy US 364	Sun Yellow US 322	Concord US 510			

Our standard upholstery shades meet the strict fire and safety regulations and are tested to Crib 5. Colours depicted may not be an exact representation of the actual material used but are the closest possible using litho printing. Please ask for a swatch for colour matching. Two-tone colour option available at no extra cost on all chairs.



# Testimonials



"I have to say that the seats have scored a hit, we have had enquiries from other theatres re the supplier as they have made such a good impression. Colour coding to each theatre eliminated any arguments about which theatre the chair belongs in, the vast colour choice was excellent, the difficulty there was selecting the colours as there was such a variety. The different specification regarding the castors and glides was really useful.

The seats are extremely comfortable and the stool tops are a generous size. They have proved to be an excellent purchase and we will certainly come back to you in the future, thank you, for your service and attention to detail."

**Theatre Sister, Spire Healthcare, Leicester**

"The products are excellent. There is a large range, with a fantastic colour selection. We therefore have a seat for everyone, including foot controlled ones for our hand and feet surgeons."

**Nuffield  
Orthopaedic  
Centre, Oxford**

"Trials were conducted with both Murray and medical staff. One distinguishing feature was the range of options available in colour, style and castor/brake options. The quality of the product can only be described as excellent."

**Golden Jubilee  
National Hospital,  
Scotland**

"Murray's chairs and stools are being used by ophthalmic surgeons and nurses in theatre. The seats are very comfortable and give good posture for our backs. The surgeons find them very supportive, as they sit throughout the surgery."

**St. Thomas' Hospital,  
London**

"My staff have found the seats extremely adjustable, very comfortable, very stable, and that the durability, wipe clean surface and good height adjustment made them a good choice for a laboratory environment."

**Wendy Cottee,  
Haematology  
Department,  
Worthing &  
Southlands  
NHS Trust**

"We tried out the full range and made the decision to purchase a selection of them for different surgeons' needs. The stools themselves are of the highest quality and are very robust. In short, my staff and surgeons are very happy with my purchase."

**Worthing Hospital,  
Sussex**

"The Great Western Hospital chose to use Murray Equipment for its theatres' surgical seating needs and have been so satisfied with the products and service that three years later have used Murray Equipment again when we opened an extra five operating suites in the new Brunel Treatment Centre.

The quality of product, alongside the after sales service given by Murray Equipment, has been exceptional and we at the Great Western Hospital would certainly recommend their services."

**Stewart Thompson, Deputy Trust Equipment Manager, Great Western Hospital, Swindon**





*Murray's*

Medical Seating

Murray Equipment Co Ltd  
Unit C1 • Ford Airfield Industrial Estate  
Arundel • West Sussex • BN18 0HY  
Tel: 01243 811881 • Fax: 01243 811855  
sales@murrayequipment.co.uk  
www.murrayequipment.co.uk

Hours of business Monday to Friday 8.30am to 4.30pm