




Committed to excellence

SL430

The **DARAY**° SL430 is an invaluable and cost effective tool. It has been specifically designed for minor procedures under local anaesthetic in health centres, GP surgeries and outpatient departments.

Ideal for minor surgery, maxillofacial, obstetrics, diagnostics and gynaecology. Also suitable for use within casualty departments and intensive care units.

The Light is designed to provide consistent, accurate light where it's needed, along with great shadow reduction and accurate colour rendering.





Features

- The use of Multiple LED's provides a high, and even light intensity throughout the entire viewing area to reduce eye strain
- A concise and self-explanatory 3 button control system on the lighthead itself
- Balanced arms for one touch adjustment and a great range of movement.
- Ideal for minor surgery and gynaecological procedures
- Removable, autoclavable handle

Environmentally friendly features

- Costs up to 75% less to run than halogen
- 50,000 hours of life from the LED light engine
- No more bulbs to change reduces waste by 95% over the light's lifetime
- Metal arm system, hardwearing and can be fully recycled
- No halogen, mercury, cadmium, lead or CFCs. No chemicals harmful to the environment.






Committed to excellence



WWW.daray.com
Tel: 0800 804 8384 Email: info@daray.co.uk
Marquis Drive, Moira, Swadlincote DE12 6EJ



**Dimension suggestion. For specific applications or environments, then mounting bracket should be installed at a height so that, when fitted, the lighthead is as near to the required position as possible

Products Available

SL430LM	DARAY SL430 LED mobile mounted minor surgery light
SL430LC	DARAY SL430 LED ceiling mounted minor surgery light
SL430LW	DARAY SL430 LED wall mounted minor surgery light
Bespoke solutions available, please call	

ABOUT DARAY

DARAY have been established in the medical, dental and veterinary markets since 1968 building exceptional and durable products.

We shall continue to innovate; aiming to excel in our endeavours to deliver reliable solutions.













Committed to excellence

www.daray.com Tel: 0800 804 8384 Email: info@daray.co.uk Marquis Drive, Moira, Swadlincote DE12 6EJ