



Code : **CH7AUTOFIT** • Reference : **3734 03** • **Lymphology**



Mobiderm Autofit - Below-knee

Size



Size	Ankle circumference	Calf circumference (to the widest part)	Circumference below the knee (popliteal fossa)
1	16 cm - 24 cm	24 cm - 35 cm	25 cm - 36 cm
2	20 cm - 28 cm	32 cm - 43 cm	32 cm - 43 cm
3	24 cm - 32 cm	40 cm - 51 cm	39 cm - 50 cm
4	28 cm - 36 cm	48 cm - 59 cm	46 cm - 57 cm



Height



Height	Height Floor / Knee (to the widest part)
Short	38 cm - 40 cm
Normal	40 cm - 42 cm
Long	42 cm - 44 cm
Very long	44 cm - 46 cm

Key Feature



Contra-indications

Do not use in contact on wounded skin, preseneting infection or acute inflammation, advanced diabted micro-angiopathy, stage 3 or 4 limb arteritis and decompensated heart faileure.

Articles

Code	Article	EAN 13	ACL
37340120100231	CHAUSS.MOB.AUTOFIT T1 COURT	34010 6014848 5	6014848
37340120100232	CHAUSS.MOB.AUTOFIT T1 NORMAL	34010 6014852 2	6014852
37340120100233	CHAUSS.MOB.AUTOFIT T1 LONG	34010 6014856 0	6014856
37340120100234	CHAUS.MOB.AUTOFIT T1 TRES LONG	34010 6014861 4	6014861
37340120200231	CHAUSS.MOB.AUTOFIT T2 COURT	34010 6014849 2	6014849
37340120200232	CHAUSS.MOB.AUTOFIT T2 NORMAL	34010 6014853 9	6014853
37340120200233	CHAUSS.MOB.AUTOFIT T2 LONG	34010 6014857 7	6014857



Code	Article	EAN 13	ACL
37340120200234	CHAUS.MOB.AUTOFIT T2 TRES LONG	34010 6014862 1	6014862
37340120300231	CHAUSS.MOB.AUTOFIT T3 COURT	34010 6014850 8	6014850
37340120300232	CHAUSS.MOB.AUTOFIT T3 NORMAL	34010 6014854 6	6014854
37340120300233	CHAUSS.MOB.AUTOFIT T3 LONG	34010 6014858 4	6014858
37340120300234	CHAUS.MOB.AUTOFIT T3 TRES LONG	34010 6014863 8	6014863
37340120400231	CHAUSS.MOB.AUTOFIT T4 COURT	34010 6014851 5	6014851
37340120400232	CHAUSS.MOB.AUTOFIT T4 NORMAL	34010 6014855 3	6014855
37340120400233	CHAUSS.MOB.AUTOFIT T4 LONG	34010 6014860 7	6014860
37340120400234	CHAUS.MOB.AUTOFIT T4 TRES LONG	34010 6014864 5	6014864
37340220100231	CHAUSS.MOB.AUTOFIT T1 COURT	31117 9740155 4	
37340220100232	CHAUSS.MOB.AUTOFIT T1 NORMAL	31117 9740156 1	
37340220100233	CHAUSS.MOB.AUTOFIT T1 LONG	31117 9740157 8	
37340220100234	CHAUS.MOB.AUTOFIT T1 TRES LONG	31117 9740158 5	
37340220200231	CHAUSS.MOB.AUTOFIT T2 COURT	31117 9740159 2	
37340220200232	CHAUSS.MOB.AUTOFIT T2 NORMAL	31117 9740160 8	
37340220200233	CHAUSS.MOB.AUTOFIT T2 LONG	31117 9740161 5	
37340220200234	CHAUS.MOB.AUTOFIT T2 TRES LONG	31117 9740162 2	
37340220300231	CHAUSS.MOB.AUTOFIT T3 COURT	31117 9740163 9	
37340220300232	CHAUSS.MOB.AUTOFIT T3 NORMAL	31117 9740164 6	
37340220300233	CHAUSS.MOB.AUTOFIT T3 LONG	31117 9740165 3	
37340220300234	CHAUS.MOB.AUTOFIT T3 TRES LONG	31117 9740166 0	
37340220400231	CHAUSS.MOB.AUTOFIT T4 COURT	31117 9740167 7	
37340220400232	CHAUSS.MOB.AUTOFIT T4 NORMAL	31117 9740168 4	



Code	Article	EAN 13	ACL
37340220400233	CHAUSS.MOB.AUTOFIT T4 LONG	31117 9740169 1	
37340220400234	CHAUS.MOB.AUTOFIT T4 TRES LONG	31117 9740170 7	
37340320100231	CHAUSSETTE MOBID.AUTOFIT T1 S	31117 9005444 3	6323165
37340320100232	CHAUSSETTE MOBID.AUTOFIT T1 N	31117 9005445 0	6323166
37340320100233	CHAUSSETTE MOBID.AUTOFIT T1 L	31117 9005446 7	6323167
37340320100234	CHAUSSETTE MOBID.AUTOFIT T1 XL	31117 9005447 4	6323168
37340320200231	CHAUSSETTE MOBID.AUTOFIT T2 S	31117 9005448 1	6323169
37340320200232	CHAUSSETTE MOBID.AUTOFIT T3 N	31117 9005449 8	6323171
37340320200233	CHAUSSETTE MOBID.AUTOFIT T2 L	31117 9005450 4	6323172
37340320200234	CHAUSSETTE MOBID.AUTOFIT T2 XL	31117 9005451 1	6323173
37340320300231	CHAUSSETTE MOBID.AUTOFIT T3 S	31117 9005452 8	6323174
37340320300232	CHAUSSETTE MOBID.AUTOFIT T3 N	31117 9005453 5	6323175
37340320300233	CHAUSSETTE MOBID.AUTOFIT T3 L	31117 9005454 2	6323176
37340320300234	CHAUSSETTE MOBID.AUTOFIT T3 XL	31117 9005455 9	6323177
37340320400231	CHAUSSETTE MOBID.AUTOFIT T4 S	31117 9005456 6	6323178
37340320400232	CHAUSSETTE MOBID.AUTOFIT T4 N	31117 9005457 3	6323179
37340320400233	CHAUSSETTE MOBID.AUTOFIT T4 L	31117 9005458 0	6323180
37340320400234	CHAUSSETTE MOBID.AUTOFIT T4 XL	31117 9005459 7	6323181
37340420100231	CHAUSSETTE MOBID.AUTOFIT T1 S	31117 9005460 3	
37340420100232	CHAUSSETTE MOBID.AUTOFIT T1 N	31117 9005461 0	
37340420100233	CHAUSSETTE MOBID.AUTOFIT T1 L	31117 9005462 7	
37340420100234	CHAUSSETTE MOBID.AUTOFIT T1 XL	31117 9005463 4	
37340420200231	CHAUSSETTE MOBID.AUTOFIT T2 S	31117 9005464 1	



Code	Article	EAN 13	ACL
37340420200232	CHAUSSETTE MOBID.AUTOFIT T2 N	31117 9005465 8	
37340420200233	CHAUSSETTE MOBID.AUTOFIT T2 L	31117 9005466 5	
37340420200234	CHAUSSETTE MOBID.AUTOFIT T2 XL	31117 9005467 2	
37340420300231	CHAUSSETTE MOBID.AUTOFIT T3 S	31117 9005468 9	
37340420300232	CHAUSSETTE MOBID.AUTOFIT T3 N	31117 9005469 6	
37340420300233	CHAUSSETTE MOBID.AUTOFIT T3 L	31117 9005470 2	
37340420300234	CHAUSSETTE MOBID.AUTOFIT T3 XL	31117 9005471 9	
37340420400231	CHAUSSETTE MOBID.AUTOFIT T4 S	31117 9005472 6	
37340420400232	CHAUSSETTE MOBID.AUTOFIT T4 N	31117 9005473 3	
37340420400233	CHAUSSETTE MOBID.AUTOFIT T4 L	31117 9005474 0	
37340420400234	CHAUSSETTE MOBID.AUTOFIT T4 XL	31117 9005475 7	