



Comfortable Respiration

F1B Full Face Mask

F1B, double-cushion facial mask, is able to provide you with excellent respiratory care. Humanized design makes the therapy as comfortable as possible.



Double cushions, more stability

F1B keeps BMC`s classic double-cushion design with more comfort and stability.

Accurate forehead adjustment

Wide adjustment range and accurate locations ensure the mask is able to fit most of the facial contours.

Upgraded headgear clips

Easy lock, quick release.


Multifunctional ports

Bilateral ports make the oxygen input or pressure titration as convenient as possible.

Product Codes

F1B Full Face Mask

- 10030162 (L)
- 10030161(M)
- 10030160 (S)



Cushion

- 60360050 (L) 60360057(M) 60360055 (S)



Thin Pad

- 60360049 (L) 60360056(M) 60360051 (S)

