



**STATIONARY
OXYGEN CONCENTRATOR OLV-10 ▲**

Advanced Features

- ▲ France imported CECA molecular sieve
- ▲ High oxygen purity 93% ± 3%
- ▲ Long lifespan more than 18000 hours
- ▲ 2 years warranty

Big LCD Display
Voice Function
24/7 Use

Specification

Specification	
Model	OLV-10
Flow Rate	0-10L/min
Purity Oxygen	93%±3%
Voltage	220V ± 22V, 50Hz ± 1Hz/110v ± 15V,60Hz ± 1Hz
Output Pressure	20kPa~60kPa
Power Consumption	550W
Noise	50dB
Net Weight	20kg
Dimension	344×306×565(mm)
Standard Accessories	1× Power cord 1× Mouth-contained nebulizer (only nebulizer model) 1× Spare first filter 1× Hepa filter 2× Nasal cannula 1× User manual 1× Warranty card


Optional Function

Nebulizer: nebulizing rate ≥ 0.15ml/min

Dual flow: Two patients can use the machine at the same time.



Nasal Cannula



Hepa Filter



First Filter



Humidifier