



Philips
Performance tubing (22
mm)

Light and flexible to use
Ergonomically designed

HH1402

New Pure White Performance Tubing (22mm)

Light and flexible

Perfect choice for smaller, lightweight masks, as it is less likely to pull on the mask and disrupt the seal.

Compatibility

- Performance tubing works well with lightweight masks

Ergonomical design

- Easy to handle

PHILIPS

Performance tubing (22 mm)
Light and flexible to use Ergonomically designed

Specifications

Design

- Compatible with: most CPAP therapy devices

Specifications

- Cuff diameter: 22mm mm
- Length: 1.82 m

HH1402/00


Highlights

Compatibility



Compatible with nearly all CPAP and BiLevel systems requiring a 22 mm non-heated hose.

Ease of use



The ergonomically designed cuff is easy to grasp and remove from the mask or machine.



Issue date 2017-05-22

Version: 5.0.1

12 NC: 0000 010 32907
UPC: 6 06959 00217 3

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com