



GammaLightTherapy.com

Thank you for purchasing a **Gamma 40 Hz Light Therapy Kit**. Your kit is *plug-and-play*, requiring no special setup. Please read the following instructions and disclaimer before use.

• • •

INSTRUCTIONS

1. Plug the 5V DC Power Adapter provided with your light therapy kit into the Gamma Light Therapy adapter.
2. Unplug to turn the lights off.
3. Enjoy your 40 Hz Gamma Light Therapy Kit.

SUGGESTIONS

- Do not use any of our products if you have photosensitive epilepsy, have a history of seizures, or experience any level of discomfort.
- We do not recommend using your light kit for more than one hour a day.
- A standard wall plug-in timer can be used to turn the lights on and off.
- Use the lights on a wall above or around a sitting room.
- Wrap the lights around an artificial plant, potted bush or tree.
- Turn the lights off at dusk so they do not interfere with sleep.

DISCLAIMER: GAMMA LIGHT THERAPY LLC. AND ITS PRODUCT MAKERS AND WEBSITE OWNERS MAKE NO GUARANTEES OR WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO PRODUCT FUNCTIONS, WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR OF NON-INFRINGEMENT OF THIRD PARTY RIGHTS. USE OF OUR PRODUCTS IS COMPLETELY AT THE USER'S RISK. OUR PRODUCTS AND STATEMENTS HAVE NOT BEEN EVALUATED BY THE FOOD AND DRUG ADMINISTRATION. THIS PRODUCT IS NOT INTENDED TO DIAGNOSE, TREAT, CURE OR PREVENT ANY DISEASE OR ILLNESS. WE DO NOT ASSUME ANY RESPONSIBILITY FOR ANY POTENTIAL ERRORS CONTAINED ON OUR WEBSITE OR PRODUCT INFORMATION PROVIDED OR ANY DAMAGES OR LOSSES RESULTING FROM THE USE OF ANY OF OUR PRODUCTS OR RELYING ON ANY INFORMATION WE'VE PROVIDED. NOT INTENDED FOR USE BY PERSONS UNDER THE AGE OF 18.