

## Home Care Medical

# Ramps

## *Threshold, Singlefold and Multifold*

Since 1974, Home Care Medical has been committed to providing the highest quality home care products and services throughout Southeastern Wisconsin. As an organization with a history of quality, responsiveness and creativity, Home Care Medical is always looking for new and innovative ways to meet the needs of our clients. That's why, Home Care Medical is pleased to offer Threshold, Singlefold and Multifold Ramps. Lightweight, easy-to-handle and set-up, our Ramps feature an anti-slip, high traction surface – perfect for scooters, wheelchairs or power chairs.

Prior to ordering a Ramp, we will be asking the following:

1. Where will you be using the ramp? (home, vehicle, etc.)
2. Will the equipment be occupied while using the ramp, and if so . . .
  - Will someone will be assisting you (and you are in a scooter, wheelchair or power chair)?
  - Or will you be moving yourself up and down the ramp on your own?
3. Where do you enter your home? (front door, side door, back door, garage, etc.)
4. How many stairs do you have to get into your home, outside and inside?
  - Are there any stairs in your yard? If yes, how many?
  - Are there any landings between the stairs?
  - If there is a landing, how long is it?
5. Would it be possible for you to measure the rise of the steps?
6. Would it be possible to measure the width of your doorway?
7. If the ramp was placed on the stairs, from the bottom step would there be enough room to have the ramp extend out for proper fit? Are there any obstacles that are close to the bottom step (bushes, flower bed, curve in sidewalk, fence, etc.)?
8. Do you have a hard time getting over the threshold of your doorway with your current mobility device? Would it be possible for you to measure the rise of the threshold?

Our Home Care Medical Personal Care Consultants are happy to assist you in choosing the perfect Ramp for your home, business or travel needs. Visit our Milwaukee, Sheboygan or West Bend Retail Stores today or call 262.786.9870 with any questions you may have. *Home Care Medical . . . helping you get on with life.*



### **Milwaukee Retail Store**

4818 South 76th Street  
Ph 414.423.8800

### **Sheboygan Retail Store**

2922 South Business Drive  
Ph 920.458.5768

### **West Bend Retail Store**

1709 South 18th Avenue  
Ph 262.957.5501

### **Corporate Headquarters**

Ph 262.786.9870  
[www.homecaremedical.com](http://www.homecaremedical.com)



## Threshold Ramps



Ramp	Rise	Length	Width	Weight Capacity	Price
EZATAEP08	Up to 1.312"	8"	32"	700 lbs.	\$77.18
EZATAEP10	Up to 1.625"	10"	32"	700 lbs.	\$82.74
EZATAEP12	Up to 2"	12"	32"	700 lbs.	\$88.10

- ✓ Accommodates thresholds up to 2" high for the perfect fit
- ✓ Designed for doorways
- ✓ Made of high-strength aluminum that will not rust or corrode
- ✓ Single-piece, formed plate with an applied, slip-resistant surface

## Singlefold Ramps



Ramp	Length	Width	Weight	Weight Capacity	Price
EZASF05	5 feet	29.5"	12 lbs.	800 lbs.	\$146.90
EZASF06	6 feet	29.5"	16 lbs.	800 lbs.	\$188.90
EZASF07	7 feet	29.5"	20 lbs.	800 lbs.	\$251.90
EZASF08	8 feet	29.5"	26 lbs.	800 lbs.	\$314.90

- ✓ Separates into two pieces for easy carrying
- ✓ Folds in half, carries like a suitcase with ergonomically designed handles
- ✓ Durable, welded construction
- ✓ Self-adjusting transition plates adjust for easy conversion from ramp to ground
- ✓ Anti-slip, high traction surface
- ✓ Full platform provides excellent stability
- ✓ Accommodates wheelchairs and scooters with various wheel configurations

## Trifold Ramps



Ramp	Length	Width	Weight	Weight Capacity	Price
EZATFAS05	5 feet	29.5"	32 lbs.	800 lbs.	\$390.50
EZATFAS06	6 feet	29.5"	38 lbs.	800 lbs.	\$457.70
EZATFAS07	7 feet	29.5"	44 lbs.	800 lbs.	\$533.30
EZATFAS08	8 feet	29.5"	50 lbs.	800 lbs.	\$690.80
EZATFAS10	10 feet	29.5"	60 lbs.	800 lbs.	\$953.30

- ✓ Separates into two pieces for easy carrying
- ✓ Folds and carries like a suitcase
- ✓ Patented hinge provides smooth operation and reduced pinch points
- ✓ Lightweight, easy-to-handle and set-up
- ✓ Made from aircraft-quality aluminum
- ✓ Self-adjusting transition plates adjust for easy conversion from ramp to ground
- ✓ Full platform provides excellent stability
- ✓ Anti-slip, high traction surface
- ✓ Accommodates wheelchairs and scooters with various wheel configurations