

DISPOSABLE BIOLOGICAL TEST PACK—EO GAS

With Instant Readout Integrator

Catalog No. BI 102

- **TEST PACK PARALLELS PERFORMANCE OF AAMI TEST PACKS**
Disposable test pack provides reliable, consistent results test after test at the lowest cost possible.
- **CONTAINS INTEGRATOR CARD**
Integrator card inserted in the test pack will not complete a color change until after biologicals have been killed. Card provides a permanent record of the success of the cycle.
- **CONTAINS SELF-CONTAINED BIOLOGICAL INDICATOR**
- **SMALL, COMPACT EASY TO USE**
The SteriTec Biological Test Pack for EO Gas sterilization with Instant Readout Integrator is assembled in a 3" x 9" pouch. Easy to store and easy to use.



ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
BI 102	30 PACKS/CASE PLUS 5 CONTROLS

SteriTec

PRODUCTS, INC.

599 TOPEKA WAY • SUITE 400
CASTLE ROCK, CO 80109 • USA
(303) 660-4201 • Fax: (303) 660-4213
Website: <http://www.steritecproducts.com>
E-Mail: info@steritecproducts.com