

Steam Biological Test Pack with Instant Readout Integrator

BI 103/INT

The Steam Biological Test Pack with Instant Readout Integrator is a steam sterilization monitor designed specifically for biological testing of pre-vacuum steam sterilizers operating at 134°C for 3 minutes or longer. The Instant Readout Integrator Card gives the operator a prediction of the biological test outcome.



Key Features

- 3 minutes or longer at 134°C
- Small, compact and easy to use
- Lead and toxic heavy metal free

Performance Characteristics

When exposed to steam sterilization cycles, the word "PASS" printed on the Integrator Card will change from purple to green, indicating an acceptable sterilization process. The biological indicator, when incubated, will indicate either growth or no growth.

Physical Characteristics

- Dimensions Card: L 101mm x W 63mm

Color Change

- Instant Readout Integrator - Purple to green
- BI Vial Process Indicator - Blue to black
- BI Growth Medium - Purple to yellow with growth

Incubation Time

24 hours at 55°-60°C

Storage

The storage conditions are printed on the product packaging. Do not store near sterilants or other chemicals.

Shelf Life

Up to 24 months. The expiry date is printed on the product packaging

Standards Classification

- Instant Readout Integrator Card: FDA Integrator
- Biological Indicator: ANSI/AAMI/ ISO 11138
- Test Pack: ANSI/AAMI ST79

SteriTec Part No.

BI 103/INT, International Steam Biological Test Pack

Packaging

30 test packs + 5 controls

FTP-044-1 Rev. G