



Product Technical Data Sheet

Form No.: FTP-056 Rev.: F

Product Name: SteriTec Cross-Checks F

Product Type: Formaldehyde Indicator

SteriTec Reorder Number: CI 118

Description: The SteriTec Cross-Checks F was designed for use in Low Temperature Steam/Formaldehyde sterilizers. The indicator check mark will change from purple to green to give visual confirmation of exposure to sterilization conditions.

Performance Characteristics:

The indicator check mark will change color when exposed to steam/formaldehyde sterilization conditions.

Color Change: Purple to Green

Indicator Stability: Excellent before and after sterilization exposure

Physical Characteristics:

Background Print Color: Text = Black, Graphic = Green (Light 7489U)

Dimensions: 3.8125" X .625"

Substrate Material: Polystyrene

Units / Bag: 250

Bags/Case: 24

Selling Unit: Case

Environmental Limitations:

Storage: Room Temperature 10° – 38°C (50° - 100°F) and 10 - 90% relative humidity

Shelf Life: The expiry date is printed on the product packaging.

SteriTec FDA 510(k) Release to Market Date: Not Released. Not for USA distribution.

Standards Classifications:

ISO 11140-1