

INTERPRETATION GUIDE

Cross-Check S - CI 134

UNEXPOSED

CROSS-CHECKS S STERILIZATION MONITOR - STEAM ISO 11140-1 CLASS 4 LOT 070401	BLUE/BLU/AZUL = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	
--	---	---

CORRECT COLOR CHANGE

CROSS-CHECKS S STERILIZATION MONITOR - STEAM ISO 11140-1 CLASS 4 LOT 070401	BLUE/BLU/AZUL = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	
--	---	---

FAIL

CROSS-CHECKS S STERILIZATION MONITOR - STEAM ISO 11140-1 CLASS 4 LOT 070401	BLUE/BLU/AZUL = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	
--	---	---



74 Inverness Drive East - Englewood, Colorado 80112

Phone: 303-660-4201 - Fax: 303-660-4213 - E-mail: info@steritecproducts.com -

Website: www.steritecproducts.com

MK-064 Rev. A