

# Cross-Checks P

## CI 115

The Cross-Checks P gives visible confirmation that the goods have been exposed to hydrogen peroxide.



### Key Features

- Process indicator for monitoring H<sub>2</sub>O<sub>2</sub> sterilization cycles
- Indicator can be used as a permanent record

### Performance Characteristics

When exposed to H<sub>2</sub>O<sub>2</sub> sterilization process, the Plasma Indicator will change from purple to a shade of pink.

### Physical Characteristics

Dimensions: L 97mm x W 14mm

### Color Change

Purple to pink

### Storage

The storage conditions are printed on the product packaging. Store in resealable bag provided away from oxidizing agents

### Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

### SteriTec Part No.

CI 115, Cross-Checks P

### Packaging

24 bags of 250