


Distributed by:
ZEPF SURGICAL
800-552-3441



Stille needle holders





Few surgical instruments are subject to more abuse than needle holders. Hence we have put great effort into producing needle holders of the highest quality. With a Stille needle holder, the needle is held firmly and safely. Most Stille needle holders are equipped with Tungsten Carbide inserts for extra durability and secure handling. Several versions are available to suit all modern suture needles. Stille needle holders offer the unique Stille feeling with a perfectly balanced and precisely aligned instrument that will work through many years of use.

Larger suture needles than the recommended sizes can damage the needle holder and invalidate the warranty. For a recommendation of needle holders according to various suture sizes, please see the “Guide for choice of needle holder in relation to suture size”.

How to test Stille Needle Holders:

First, clamp the needle in the jaws of the needle holder, then lock the instrument in the second ratchet tooth.

If the needle can be turned easily by hand, the instrument should be set aside.

When the instrument is new, it should hold securely on the first ratchet tooth for a considerable time.

While holding the rings of the instrument in each hand, open the instrument as far as possible. If the box lock has considerable play and it is very loose, this will cause misalignment. If this occurs, the instrument is faulty and should be repaired.

Cleaning and sterilizing of a needle holder should always be carried out with the ratchet open.



Guide for choice of needle holder in relation to suture size

| Suture size | 1/0 | 2/0 | 3/0 | 4/0 | 5/0 | 6/0 | 7/0 | 8/0 | 9/0 | 10/0 | 11/0 |
|---------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|------|
| 104-15542 | | | | | x | x | | | | | |
| 104-15545 | | | | | | | x | x | | | |
| 104-15546 | | | | | | | x | x | | | |
| 104-15550 | | | | x | x | x | | | | | |
| 104-15551 | | | | x | x | x | | | | | |
| 104-15553-17 | | | | | x | x | x | | | | |
| 104-15556 | | | | | x | x | | | | | |
| 104-15557 | | | | | x | x | | | | | |
| 104-15558 | | | | | x | x | | | | | |
| 104-15559 | | | | | x | x | | | | | |
| 104-15560 | | | | x | x | x | | | | | |
| 104-15560-2 | | | | x | x | x | | | | | |
| 104-15561-141 | x | x | x | x | | | | | | | |
| -161 | x | x | x | x | | | | | | | |
| -181 | x | x | x | x | | | | | | | |
| -201 | x | x | x | x | | | | | | | |
| -251 | x | x | x | x | | | | | | | |
| 104-15565-151 | | | | x | x | x | | | | | |
| -181 | | | | x | x | x | | | | | |
| -201 | | | | x | x | x | | | | | |
| -231 | | | | x | x | x | | | | | |
| 104-15566-151 | | | | x | x | x | | | | | |
| -181 | | | | x | x | x | | | | | |
| -201 | | | | x | x | x | | | | | |
| -231 | | | | x | x | x | | | | | |
| 104-15569 | | | | x | x | x | | | | | |
| 104-15571 | | | x | x | x | x | | | | | |
| 104-15573-161 | | | | x | x | x | | | | | |
| -181 | | | | x | x | x | | | | | |
| -201 | | | | x | x | x | | | | | |
| 104-15585-161 | | | | | | x | x | x | | | |
| -181 | | | | | | x | x | x | | | |
| -221 | | | | | | x | x | x | | | |
| -251 | | | | | | x | x | x | | | |
| 104-15586 | | | | | | x | x | x | | | |
| 104-15586-201 | | | | | | x | x | x | | | |
| | | | | | | | | | | | |
| | | | | | | | | | | | |
| 104-15621-179 | x | x | x | | | | | | | | |
| -209 | x | x | x | | | | | | | | |
| 104-15777 | | | | | | x | x | x | | | |
| 104-15778 | | | | | | x | x | x | | | |

Distributed by:

ZEPF SURGICAL

800-552-3441



Distributed by:
ZEPF SURGICAL
 800-552-3441



RAPP-STILLE

A 13 cm knot-holding forceps with smooth jaws.

* 104-15140- 13 cm 5"



DERF

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 5/0 to 6/0 suture.



* 104-15542- TC jaws 12 cm 4.75"



GILLIES

Needle holder with scissors. Ergonomic design on shanks and finger rings. One rather large ring handle is set at an angle for fitting the base of the thumb. Without ratchets. With or without TC jaws. The scissor part enables suturing and cutting with one instrument. Intended for 4/0 to 6/0 suture.

* 104-15550- With Scissors 16 cm 6.25"
 104-15551- With Scissors TC jaws 16 cm 6.25"




* Photo shows this instrument

NEEDLE HOLDERS

WEBSTER DELICATE

Very delicate needle holder for very fine sutures. Smooth jaws or TC jaws. Intended for 7/0 to 8/0 suture.



- * 104-15545- Ultra delicate TC-sprayed jaws 12 cm 4.75"
- 104-15546- Ultra delicate Smooth jaws 12 cm 4.75"



SENNING

Delicate needle holder with small, short-cupped cross-serrated jaws for fine sutures. Shanks flared out for comfortable grip. Ideal for 5/0 to 7/0 suture.


- * 104-15553-17 Delicate 17 cm 6.75"



WEBSTER/KLEINERT

Delicate needle holder for fine needles. With smooth jaws or TC jaws. Intended for 5/0 to 6/0 suture. The Webster-Kleinert (104-15557) has offset bows for a comfortable grip.

- 104-15556- Delicate Smooth 12 cm 4.75"
- * 104-15557- Delicate Smooth 12 cm 4.75"
Offset bows
- 104-15559- Delicate Serrated TC jaws 12 cm 4.75"





* Photo shows this instrument

WEBSTER-STILLE

Delicate needle holder for fine needles. With smooth TC jaws. Intended for 5/0 to 6/0 suture.


- * 104-15558- Smooth TC jaws 12 cm 4.75"



HALSEY

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb. For slightly larger needles than the Webster needle holder. Intended for 4/0 to 6/0 suture.


- * 104-15560- TC jaws 12.5 cm 5"
- 104-15560-2 Ergo TC jaws 12.5 cm 5"



MAYO-HEGAR

Commonly used general purpose needle holder with TC jaws. Intended for 1/0 to 4/0 suture.

- * 104-15561-141 TC jaws 14 cm 5.5"
- 104-15561-161 TC jaws 16 cm 6.25"
- 104-15561-181 TC jaws 18 cm 7"
- 104-15561-201 TC jaws 20 cm 8"
- 104-15561-251 TC jaws 25 cm 10"






* Photo shows this instrument

CRILE-WOOD

General purpose needle holder with TC jaws and two different jaw widths. Ideal for 4/0 to 6/0 suture.


| | | | |
|-----------------|----------|---------|----------|
| 104-15565-151 | | TC jaws | 15 cm 6" |
| 104-15565-181 | | TC jaws | 18 cm 7" |
| 104-15565-201 | | TC jaws | 20 cm 8" |
| 104-15565-231 | | TC jaws | 23 cm 9" |
| * 104-15566-151 | Delicate | TC jaws | 15 cm 6" |
| 104-15566-181 | Delicate | TC jaws | 18 cm 7" |
| 104-15566-201 | Delicate | TC jaws | 20 cm 8" |
| 104-15566-231 | Delicate | TC jaws | 23 cm 9" |



BROWN

Needle holder with fine serrated, TC-sprayed rounded jaws, intended not to harm suture. Suitable for 4/0 to 6/0 suture.



| | | |
|--------------|-----------------------------|------------|
| * 104-15569- | TC sprayed, Diamond jaws | 14 cm 5.5" |
|--------------|-----------------------------|------------|



NEW ORLEANS

Short tapered TC jaws, one jaw straight serrated and the other cross-serrated. Intended for 3/0 to 6/0 suture.

| | | |
|--------------|-------------------------------------|----------|
| * 104-15571- | TC jaws, serrated/cross serrated | 18 cm 7" |
|--------------|-------------------------------------|----------|





* Photo shows this instrument

RYDER

Needle holder with very narrow, non-tapered, serrated TC jaws. Popular in vascular surgery. Intended for 4/0 to 6/0 suture.



- * 104-15573-161 TC jaws 16 cm 6.25"
- 104-15573-181 TC jaws 18 cm 7"
- 104-15573-201 TC jaws 20 cm 8"



STILLE NEEDLE HOLDER

TC-sprayed delicate jaws. For fine suture 6/0 to 8/0.



- * 104-15585-161 TC sprayed, Diamond jaws 16 cm 6.25"
- 104-15585-181 TC sprayed, Diamond jaws 18 cm 7"
- 104-15585-221 TC sprayed, Diamond jaws 22 cm 8.75"
- 104-15585-251 TC sprayed, Diamond jaws 25 cm 10"



STILLE NEEDLE HOLDER

Delicate needle holder with TC jaws. Ideal for 6/0 to 8/0 suture.

- * 104-15586- TC jaws 18 cm 7"
- 104-15586-201 TC jaws 20 cm 8"





* Photo shows this instrument

MATHIEU

Needle holder with straight TC jaws. Plier type handle with ratchets that allows the surgeon to open and close the needle holder rapidly. Slightly longer jaws than Metzenbaum needle holder. Intended for 1/0 to 3/0 suture.


- * 104-15621-179 TC jaws 17 cm 6.75"
- 104-15621-209 TC jaws 18 cm 7"



KAYE

Delicate needle holder with short, fine TC jaws and sharp scissors. The scissor part enables suturing and cutting with one instrument. Ideal for 6/0 to 8/0 suture.

- * 104-15777- Delicate, Serrated 12 cm 4.75"
 with scissors TC jaws
- 104-15778- Delicate, Smooth 12 cm 4.75"
 with scissors TC jaws



* Photo shows this instrument