



PLATINUM LID CHAIN Implant




(enlarged to show detail & range of motion)


A flexible chain implant made of platinum adapts better to the changing radius of the tarsus during movement of the globe.

Platinum has a higher density than gold, so a smaller implant of the same weight can be utilized.

Available in eight sizes from 0.6gm to 2.0gm.



PLATINUM CHAIN IMPLANT











STANDARD GOLD IMPLANT



THE BERGHAUS/SCHROM PLATINUM LID-CHAIN is an implant used to correct lagophthalmus. The implant is in the form of a chain with several movable links. They are sterile packed, and available in eight sizes.

PLATINUM LID CHAIN Implant

.6 gm	28-1506*	
.8 gm	28-1508*	
1.0 gm	28-1510*	
1.2 gm	28-1512*	
1.4 gm	28-1514*	
1.6 gm	28-1516*	
1.8 gm	28-1518*	
2.0 gm	28-1520*	

TESTWEIGHTS FOR LID-CHAIN
(0.6-2.0 gm) 28-1550

ADHESIVE STRIPS FOR
TESTWEIGHTS 28-1580
(box of 50)

* Indicates Sterile Item

“ A certain percentage of undesired effects and complications are known to occur when upper-lid solid gold implants are used in facial palsy. The flexible platinum chain implant adapts better to the changing radius of the tarsus with the movement of the globe and reduces the complication rate. These platinum chains are smaller due to their higher density.”

Arch Facial Plast Surg. 2003;5:166-170