#### **EVEN EASIER TO LOAD**

We redesigned our best-selling model and broke it down into 3 easy to load pieces, all capable of fitting in the back of any vehicle No need for specialty transport vehicles, vans, or expensive lifts.





## **Performance Specifications**

**Driving Range** Maximum Speed Drive System **Motor Type** Incline Capability

15.5 mi / 25 km 4.5 mph / 7 kph Two Motor Drive System 250 Watt Brushless 12.0° / 21.3%

Turning Radius Stopping Distance Obstacle Clearance **Braking System** Max Weight Capacity

5 ft (at full speed)

Electromagnetic

### **Dimensional Specifications**

**Unfolded Dimensions** Chair Gross Weight Chair Net Weight **Seating Width** Seat Cushion Size

44" L x 23" W x 40" H 65 lbs 58 lbs (no batteries) 19.5" / 23" 17" W x 18" D x 3" H

Folded Dimensions Front Wheel Diameter Rear Wheel Diameter Floor to Seat Height **Back Support Height** 

33" L x 23" W x 13" H 12" 19.5" 20.5"


## Other Specifications

Battery Type **Batteries Provided Battery Watt Hour Battery Charge Cycles** Tires

Lithium-ion (24V 6Ah) Two (2) 144 Wh (Airline Approved) ~1,000 Non-Puncture

Charger Input Lap Belt

5 to 7 hours 100-240V, 50/60HZ(international) 29.4V 2.0A Provided Built-In





CITY 2 PLUS Foldable Power Chair



# **Live Freely with Travel Buggy Electric Wheelchairs**

Convenient to use. Quick to fold. Easy to transport.





travelbuggy.com

contact@travelbuggy.com



# Your New Best Friend for Personal Mobility

You've waited long enough. It is time to see your friends and family, explore new places, maybe even travel to far off destinations, and it is all possible in our new CITY 2 PLUS. We've made our Best Selling CITY model even more convenient for you or your loved ones to fold, load, and go. It is time to get out there again. It is time to join the Travel Buggy family.



# **Easy to Fold and Transport**

Folds quickly like a lawn chair to the size of a suitcase for convenience transport in any vehicle or for storage.



1-Step Fold Sequence



Net weight 58 lbs



Fits any vehicle



## **Will Travel!**

Travel Buggy electric
wheelchairs do not require any
additional travel batteries or
approvals like other "Travel
Chairs". Our LiPo batteries are
rated for safe airline travel by
TCCA (CA), FAA (US) and
various other aviation
authorities.





## **Intelligently Designed**

Our Travel Buggy electric wheelchairs are intuitive to drive and convenient to manage.



Multi-position reclining back rest (up to 135°) to relieve back pain & ensure comfortable all-day use



Batteries are easy to remove for airline travel



Entry and exit is made easier by a fold-away footrest & raiseable arms for optional side entry



Do more indoors and outdoors with our CITY 2 PLUS model. Designed to handle what the outdoors can bring plus works great in tight spaces.



Use Indoors & Outdoors



Top Speed



Long Driving Range



## Make Your Travel Buggy Unique To You

We offer more accessories than anyone else. From essentials like phone holders & extra storage to rehab accessories like head rests and leg extensions. We can accessorize your CITY 2 PLUS power chair to your desire.

















The information contained herein is correct as at the time of printing, we reserve the right to alter specifications without prior notice