



THERMOSTATIC DEVICES

Dry incubator for antibiotics testing

www.milkotester.com
 www.milkanalyzer.eu
 email: office@milkoteter.com
 tel.: +359 888 259 646

Prevents the risk to receive milk



in which milk acidic bacteria is developing



SPECIFICATION

Model name	Compatible tests
Milkotester Rosa-W2	Charm® 
Milkotester CH-W12	Charm® 
Milkotester DC-W14	Delvotest®  Copan® 
Milkotester ET-W24	Eclipse® 
Milkotester BS-W8	Beta Star® 



- Wells
- LED Display
- START/SET button
- DOWN/STOP button
- UP button
- ON/OFF buton
- Indicator lights

High accuracy,
 easy to use and
 stable
 temperature



Stainless steel cover box,
 aluminium plate with different
 thickness (12-30mm) and
 different number of wells in
 dependence of the model

The devices can be connected
 on the milk transportation
 tank vehicle





THERMOSTATIC DEVICES

Dry incubator for antibiotics testing

www.milkotester.com
www.milkanalyzer.eu

Prevents the risk to receive milk



in which milk acidic bacteria is developing

Key features

Time range:
0 to 244sec.
0 to 600min.

Temperature range:
+30°C to +80°C

Temperature accuracy:
(± 0.5°C)

Two factory set working models

Low power consumption

Compact design

Simple in maintenance, setting and installation



ET - W24

BS - W8



Rosa - W2



DC - W14

CH - W12