

Internal printer (option) can print instantly the results from the milk analyzer or create a paper report for a client by number and liters




pH measurement (option) is used to determine the acidic or caustic composition of the milk sample



Conductivity meter (option) is used for determination of salts in the milk sample.

Accumulator (option) can supply electricity to the milk analyzer in field work up to 8 hours without charging.



Flash drive (option) is used for collecting and transferring measurement data from the milk analyzer.



Second pump (option) is used for auto cleaning and rinsing the milk analyzer with cleaning solution.

