



Product
Catalog
2019



 **MDC**[®]
DENTAL

Our Company



In dentistry and all of its branches there exists a science involving technology, art, skills and the great desire provide function and an aesthetically pleasing result to a patient.

At MDC, our odontological advances in manufacturing are backed up by three very important factors. Experience, Knowledge and Servant hood. These three factors combine to form our core philosophy of providing excellent customer service to our clients along with the best products on the market. We are committed to the trust that you give us by choosing us over anyone else.

Our company has over 150 employees that make it possible to be one of our nation's leading companies in the dental industry. Our research laboratories have one-of-a kind scientific equipment, the best accredited & specialized chemical & scientific engineers working to bettering our products each and every day.

Our partnerships with well-known national and international institutes and companies, allow us to support future project as well as give us total and permanent control of all production procedures and testing.

There are dental supply companies, and then there is MDC Dental – We welcome to opportunity to show you our core philosophy at work.

The MDC Dental Team



Lab Products

Acrylics	3
Waxes	13
Teeth	16
Gypsums	20
Separating film	22
Investments	21
Reversible hydrocolloid	22
Pastes and polishings	23
Articulators	24

Clinic Products

Alginates	25
Dental dam	30
Dental enamel cleaning	31
Dental adhesive	32

Clinical products

Adhesive paste for dental prosthesis	32
Alginates	25
Cleasing effervescent tablets	33
Dental dam	30
Dental enamel cleaning kit	31
Max Print measure mix	29
Mixing bowl	29
Powder adhesive for dental prosthesis	32
Tissue conditioning	11

Laboratory products

Acrilycs	3
Articulators	24
Investment	21
Plaster	20
Polisher	23
Resin teeth 2 layers	16
Reversible Hydrocolloid	22
Separating liquids	20
Shade guide Acrilycs Nic Tone	10
Tissue conditioning	11
Waxes	13

Brand index

Antivet	31
BioJel	29
Confident	32
Kelvin	21
Magnum	20
Max Print	25
New Flex	22
NewTek	16
Nic Tone acrylics	3
Nic Tone dental dam	30
Polycril	23
Polyshine	23
Rhino articulators	24
Rhino polishing	23
Rogson Wax	13
SepyA	20

Images are for illustration purposes only



Nic Tone Acrylics



By offering a wide range of products and solutions for the clinician and lab - Nic Tone® is the most recognized brand in the market

- Pink and transparent tones
- Tooth Colors
- Self - Setting for orthodontics
- Color liquids
- High Impact
- Monomers
- Flowable resin
- Self curing resin
- Tissue Conditioner



Nic Tone
FORZZA
High impact

Clear and pink tones

An Acrylic PMMA based polymer with polymethyl methacrylate. Having superior characteristics (Cross Linked) that reduce the risk of prosthetic porosity and increase resistance to impact and abrasion, this is by far one of the best products on the market today.

Applications

- Partial and full denture bases
- Removable partial prosthetic denture bases
- Relines
- Dental implants
- Repairs (self-curing)
- Occlusal Guards (transparent)

Characteristics

A varied and reliable color scheme

Advantages

- Easy to work with
- Applicable to any polymerization technique
- No porosity, eliminates changes in color and odor
- No bacterial plaque build up



	CLEAR	R1V LIGHT PINK	R2V MEDIUM PINK	R6V MEDIUM PINK EXTRA VETA	R7V PINK DARK
COLD CURE					
Jar 3.17 oz	107-3011	107-3021	107-3031	107-3061	107-3071
Jar 1.1 lb	107-3013	107-3023	107-3033	107-3063	107-3073
Jar 5.51 lb	107-3014	107-3024	107-3034	107-3064	

HEAT CURE				
Jar 1.1 lb	105-1013	105-1023	105-1033	105-1063
Jar 5.51 lb	105-1014	105-1024	105-1034	105-1064

POUR TYPE					
Jar 1.1 lb		114-5020	114-5030	114-5070	114-5080
Jar 5.51 lb		114-4020	114-4030	114-4070	114-4080

Note: Tones just for reference

Tooth shades

Polymethyl Methacrylate PMMA

Applications

Heat-curing: In making pontic teeth and aesthetic veneers

Self-curing: In making temporary bridges and crowns

Characteristics

Vita® shades

A1 A2 A3 A3.5 A4 B1 B2 B3 C2

New Hue shades

59 61 62 65 66 67 69 81



Coldcure	A1	A2	A3	A3.5	A4	B1	B2	B3	C2
Jar 1.23 oz	107-3210	107-3220	107-3240	107-3230	107-3250	107-3260	107-3270	107-3280	107-3290
Jar 6.70 oz	107-3212	107-3222	107-3242	107-3232	107-3252	107-3262	107-3272	107-3282	107-3292
Jar 1.1 lb	107-3213	107-3223	107-3243	107-3233	107-3253	107-3263	107-3273	107-3283	107-3293

Coldcure	59	61	62	65	66	67	69	81
Jar 1.23 oz	107-3100	107-3120	107-3130	107-3140	107-3150	107-3160	107-3170	107-3180
Jar 6.70 oz	107-3102	107-3122	107-3132	107-3142	107-3152	107-3162	107-3172	107-3182
Jar 1.1 lb	107-3103	107-3123	107-3133	107-3143	107-3153	107-3163	107-3173	107-3183

Basic VITA® shades 4 x 1

Basic Kit of colors for provisional development

Applications

Temporary crowns and bridges

Characteristics

Universal Vita® shades coldcure A1 ● A2 ● A3 ● B1 ●

Advantages

- Shades similar to the final result
- Easy to store
- Easy to use
- Monomer and dropper included

Basic Vita® shades 4 x 1	106-3003
--------------------------	----------

4 jars/1.2oz each

1 bottle coldcure liquid 4.2oz

1 dropper

1 storage rack



Custom Tray Materials

A self-curing resin material for custom trays that's easy to use and manipulate. Very stable material that does not distort after curing. Also available in Green Powder. Use our Nic Tone self-curing liquid to process.

Applications

- To make individual custom trays and bases

Characteristics

- Polymer self-curing acrylic
- Mix with self-curing monomer Nic Tone
- Green color

Advantages

- Easy to handle
- Produces individual custom trays for uniform impressions
- Does not stick
- More economic than any self-curing acrylic

Jar 1.1 lb	108-1013
Jar 5.5 lb	108-1014



Orthodontic self-curing

The most beautiful orthodontic acrylic!
Formulated with the most vibrant, attractive
and complete color line available.



Characteristics:

- Ultra fine powder 120 mesh
- Non porosity acrylic



	Fluorescent yellow 1	Yellow Glitter 2	Crystal Blue 3	Opaque White 4	Blue Glitter 5	Fluorescent Orange 6	Orange Glitter 7	Lilac Glitter 8	Opaque Black 9	Turquoise Glitter 10
1.23 oz	106-2010	106-2020	106-2030	106-2040	106-2050	106-2060	106-2070	106-2080	106-2090	106-2100
6.7 oz	106-2012	106-2022	106-2032	106-2042	106-2052	106-2062	106-2072	106-2082	106-2092	106-2102
1.1 lb	106-2013	106-2023	106-2033	106-2043	106-2053	106-2063	106-2073	106-2083	106-2093	106-2103

	Crystal Red 11	Red Glitter 12	Pink Glitter 13	Fluorescent Pink 14	Pearl Pink 15	Crystal Green 16	Fluorescent Green 17	Green Glitter 18	Violet / Violet Glitter 19 / 19D	Opaque Blue 20
1.23 oz	106-2110	106-2120	106-2130	106-2140	106-2150	106-2160	106-2170	106-2180	106-2190 / 106-2200	106-2210
6.7 oz	106-2112	106-2122	106-2132	106-2142	106-2152	106-2162	106-2172	106-2182	106-2192 / 106-2202	106-2212
1.1 lb	106-2113	106-2123	106-2133	106-2143	106-2153	106-2163	106-2173	106-2183	106-2193 / 106-2203	106-2213

Pro kit Ortho Powder and Liquid

The finest orthodontic acrylic colors!
The formulation is bright and attractive
with a complete line of colors.

Applications:

- Development of Orthodontic and orthopedic retainers

Characteristics

- Extra fine polymer provides a better finish
- Twenty bright colors
- Pro kit comes with all shades and liquid

Advantages

- Easy to use extra fine acrylic
- Versatile techniques due to excellent consistency
- Rainbow glitter



Pro kit ortho	106-3010
---------------	----------

1 Fluorescent yellow	2 Yellow Glitter	3 Crystal Blue	4 Opaque White	5 Blue Glitter	6 Fluorescent Orange	7 Orange Glitter	8 Lilac Glitter	9 Opaque Black	10 Turquoise Glitter
11 Crystal Red	12 Red Glitter	13 Pink Glitter	14 Fluorescent Pink	15 Pearl Pink	16 Crystal Green	17 Fluorescent Green	18 Green Glitter	19 Violet Glitter	20 Opaque Blue

Professional Liquid Color Kit

Applications

- Fabrication of orthodontic and orthopedic appliances

Characteristics

- Color monomer included to mix with self-curing polymer

Advantages:

- 6 basic colors
- Versatility in color mixing



Liquid Color Kit	141-0010
------------------	----------

Includes:

6 colors in bottles of 4.2 oz

1 Transparent self-curing polymer of 6.7 oz

1 glitter 2.4 oz

1 dropper

Self-Curing liquid: 125 ml

Blue	141-0030
Purple	141-0040
Green	141-0070
Yellow	141-0020
Orange	141-0050
Red	141-0060

Forzza High Impact Acrylic



Superior impact resistance due to its unique innovative formula. Made for extreme use with a natural appearance.

Applications

Base for partial and full dentures. Aesthetic coating for prosthetic implants.

Characteristics

- Easy to use
- High resistance to fracture
- Pleasant pink tone
- Reliable shade

Advantages

Recommended for patients with bruxism
Extremely resistant

Forzza Kit	115-1010
Includes:	
1 Jar (polymer) 6.7 oz	
1 Bottle (monomer) 4.2 oz	
6 Bags	
1 instructive	



Tissue Tone Technique

Chromatic Acrylic system for custom prosthesis provide the best natural look. (Heatcure kit)

Application

- For the fabrication of prosthesis with a natural look

Tissue tone kit	109-1010
-----------------	----------

Contain:

5 jars 1.2 oz

3 dispensers

1 dropper

10 dosifier tips

		Soft pink		Light pink		Violet		Capilar pink		Artery red
Jar 1.2 oz	109-1050	109-1030	109-1060	109-1040	109-1020					




Complete Pedestal Shade Guide

Applications

Nic Tone Color Guide: Forzza, Custom Trays, Transparent, Orthodontic, Pink Nic Tone and Tooth shades.

Complete Pedestal Shade Guide	116-0002
-------------------------------	----------



Pink color tone shade guide

Applications

Acrylic color guide for Nic Tone Pink and Forzza

R1V	R2V	R3V	R4V	R5V	R6V	R7V	R8V	R9V	Forzza

Note: Tones just for reference

10 tones Pink shade guide	116-0001
---------------------------	----------



Inlay Pattern Resin



Nucleo Lay will radically change your casting techniques!

Applications

- Casts for the reconstruction of post/stumps.
- Incisal guides
- Calibrating and fastening

Characteristics

- Dimensional stability
- Total absence of solid waste to be incinerated
- Fast working time

Advantages

- For lab and clinical use
- Minimum contraction
- Calcinable



Nucleo Lay Kit <i>Includes:</i> 1 Resin bottle 1.4 oz 1 Monomer bottle 1.3 oz 1 separator bottle 0.5 oz 1 brush 1 dropper 10 plastic posts 1 instructive	119-0010
1 Resin bottle 1.4 oz	119-0030
1 Monomer bottle 1.3 oz	119-0020
16 plastic posts	

Tissue conditioning

NIC TONE
Softy

An advanced soft material for relines with outstanding durability.
Coldcure up to 8 months and Heatcure (Lab processing)
up to 24 months.

Applications

- Material for direct and indirect relines
- Tissue conditioning

Characteristics

- Resilient for partial or full denture liner
- Long lasting up to 24 months if heat cured
- Lasts up to 6 months if self cured

	Self curing	Heat curing
Full Kit Powder Liquid Test tube Instruction Manual	120-0020	120-0010
Replacement Kit Powder Liquid Instruction Manual	120-0040	120-0030
1 Copolymer bottle 2.8 oz	120-0104	120-0102
1 Monomer bottle 1.3 oz	120-0103	120-0101



Liquid Acrylic

Heatcure and Coldcure monomer

Applications

To be used with polymers

Characteristics

- Clear liquid

Advantages

- Glycol ether added to monomer to induce cross link polymerization
- Sold in dark glass containers and metal containers to preserve its original properties and its shelf life.

	Heatcure	Coldcure	Pourtype
1/4 gal (32 oz)	137-0050	138-0050	114-2040
Container 1 gal.	137-0070	138-0070	114-2030



*Presentation available in USA

All Purpose Base Wax

Wax for all kinds of dental applications (also baseplate). Excellent quality made with soft and flexible materials, works in extreme temperatures. The entire line has ideal characteristics for demanding laboratory jobs.

Applications

All purpose, setup and base

Characteristics

Consistency:

Medium all purpose Medium Pink

Medium Ethnic Dark Pink

Hard Pink

Soft Red

Advantages

- Easy to use
- Non sticky
- Flexible
- Holds teeth firmly
- Wax is fast and efficient
- Great for samples and trials
- Clean extraction, leaves no residue



Complies with ISO 15854

	PINK 1	DARK PINK 2	PINK 3	RED 4	PINK 5
Consistency	Medium all purpose	Medium ethnic	Hard	Soft	Extra hard
Box with 1/2 lb	148-0061	148-0064	148-0062	148-0060	148-0063
Box with 5 lb	148-0070	148-0073	148-0071	148-0074	148-0072

Anatomic U shape bite block

Its unique shape is designed to get the best results in making impressions of teeth. Anterior section is deeper than the posterior section, avoiding unnecessary carving or adjustments.

Applications

Partial and full Maxillo-mandibular impression


Characteristics

Consistency / Color:

- Soft Red
- Medium Pink
- Hard Light Pink

Advantages

- Easy to use
- Non Sticky
- Excellent flexibility
- Holds teeth firmly
- Fast and efficient wax
- Great for test samples
- Clean extraction, leaves no residue



	RED	PINK	PINK
Consistency	Soft	Medium	Hard
Box with 72 pieces	148-0103	148-0101	148-0102

Universal Palate Pattern



Flexible and adaptable anatomic wax palate

Applications

- For making Acrylic and flexible dentures

Characteristics

- All purpose wax
 - Flexible and adaptable
 - Wax palate functional and natural
 - Quick phonic adaptation
 - Quick application to patients ·
- Quick applications in the laboratory



CE

Palate Box w/50 pcs.	148-0080
Strip Vestibular Box w/100 pcs.	148-0090

Utility Wax



For post-damming and build up purposes

Applications

- Impressions
- Block retentive areas
- Brackets (clear)

Characteristics

Round strips
Red / white

Advantages

- Dimensional stability
- Distortion free
- Easy to separate
- Soft
- Excellent consistency



	RED	WHITE
Box of round strips 11.81 oz	148-0031 	148-0030 

CE

Round wax, wire spools

Easy setup and perfect hardness



Applications

- Easy to form wax wire used for producing various bars and sprues for partial dentures, crowns and bridge work.

Characteristics

- Firm and moldable

Advantages

- Variety of calibers

#6	#8	#10	#12	#14	#16
300 g / 10.58 oz	240 g / 8.46 oz	225 g / 7.93 oz	225 g / 7.93 oz	240 g / 8.46 oz	225 g / 7.93 oz
148-0051	148-0052	148-0053	148-0054	148-0055	148-0056

#18	#20
10 g / 0.35 oz	10 g / 0.35 oz
148-0058	148-0057



CE

Probase

Thermoformable rigid base

Uses:

For base elaboration.

Characteristics:

5" X 5" sheet
Gauge 0.060 (0.5 mm)


Advantages:

- Fast (use vacuum-form, process it and cut)
- Clean (Its use avoids wax work on the base.)
- Stable (For the maxillomandibular measurement, without changes and in a predictable way).

Box w/25 pc	122-0020
-------------	----------



Upper anteriors



Lower anteriors



A23

43.5




30M




A23

34.0



A24

46.0



30M

29M




A24

40.0



A25

49.0



30M

29M

32F

33F




A25

40.0



A26

52.0



32M

32L

31M

32F

33F

31Z




A26

41.5



3D

48.0



32M

32L

32F

33F

31Z




3D

39.0



3N

48.0



30L

32L

31Z



3N

38.0






Upper anteriors

Lower anteriors



50.0



32M 32L 31M 32F 33F 31Z



39.0



53.0




32L 33M 33L 32F 33F



42.0



53.5




32M 34M 34L 33M 32F 33F 31Z



43.5



57.0



34M 34L 33M 32F 33F



44.0



45.0



34M 34L 33M 32F 33F



42.0



55.0



34M 34L 33M 32F 33F




42.0



63.0








34M 34L 33M 32F 33F



44.5



Posteriors



newtek *Acrylic resin teeth*

The highest quality material and finish as well as great resistance

- Crosslinked
- Natural fluorescence
- Computerized pigment
- Long durability
- Various anterior and posterior
- Molds available in 33°, 20°, 10°



Color guide



New Hue* shades
59, 61, 62, 65, 66, 67, 81

69 and 77 ask for availability



Vita** shades
A1, A2, A3, A3.5, A4, B1, B2, B3, C2

Posteriors: Box with 20 pieces
Anteriors: Box with 30 pieces

* New Hue is a registered trademark of DENTSPLY International Inc.
** Vita is a registered trademark of VITA International.

Child teeth



Magnum plaster type III



The Quality of our materials and manufacturing processes make MDC Dental's plasters the most reliable product in the market.

Manufactured under ISO and ADA specifications & ISO 9001 - 2008 standard

Applications:

- Study models
- Antagonist bases
- Work models
- Bases

Characteristics

- Class III super stone
- Regular setting
- Extra fine powder
- Optimal hardness

Advantage

- Smooth surface
- Creamy smooth gypsum
- High detail
- Low setting expansion
- Dimensionally stable

Physical properties

Setting expansion
Compression resistance
Setting time
Water/Powder ratio

Normal

0.15%
390 kg/cm²
5'30" to 8' min.
24 cc / 100 gr

- Study models
- For bases and antagonist bases

Ortho

0.15%
320 kg/cm²
5'30" to 8' min.
36 cc / 100 gr

- Models for orthodontics and antagonists



NORMAL	OCHRE	WHITE	BLUE
Box with 22 lb	129-0013	129-0033	129-0023



ORTHO	WHITE
Box with 22 lb	129-0066

Paris type II

Controlled expansion for perfect joints

Applications

- Incased in plaster
- Study models
- Antagonist bases
- Work models
- Base

Characteristics

- Smooth surface
- Regular setting
- Extra fine powder
- White
- Dimensionally stable

Physical properties

Setting expansion
Compression resistance
Setting time
Water/Powder ratio

0.20%
230 kg/cm²
5'30" a 8'
41 cc / 100 gr



Bulk with 22 lb	129-0091
-----------------	----------

Premium investment for chrome cobalt casting



Applications

- Investment casting in chrome cobalt

Characteristics

- Liquid to control expansion
- Fine micro grain
- High Precision
- Clean casting



Bucket of 25 bags / 2.2 lbs each

134-0010

Phosphate investment for nickel-chrome alloy



Applications

- Investment wax patterns for casting in ceramic alloys

Characteristics

- Liquid to control expansion
- Controlled expansion
- High temperature resistance



Kit: 10 bags w/2.2 lbs each of powder

132-0015

Reversible Hydrocolloid

Applications

- For plaster and investment models duplications

Characteristics

- Gel in colloidal particles
- Rapid gel

Advantages

- Maximum accuracy
- Highly reusable
- All purpose (investment and gypsum)

**SUPER
ALL PURPOSE**

Bucket 20 lb	125-0021
--------------	----------

NEW FLEX



Sepya® Gypsum and acrylic separator

Applications

Liquid gel designed to prevent the plaster from sticking to the acrylic

Advantages

- Clean separation
- Easy removal of the plaster from the acrylic
- Avoids altering the precision of the thicknesses
- Flexible to 0.05 microns

Bottle 32 oz (1/4 Gal)	126-0030
Container 1 Gal	126-0040



Degreasing liquid eliminates superficial tension

Applications

- Cleans oil released from wax
- Conditions the wax patterns
- Eliminates superficial tension
- Eliminates the majority of air bubbles from wax patterns/casts

Characteristics

- Spray bottle

Advantages

- Easy to apply

Spray bottle 4.2 oz	127-0010
---------------------	----------

Finitto



Polish for acrylics, resins, and flexible materials

polycril

Applications

- Acrylics, resins, and flexible materials.

Characteristics

- 100% natural
- Produces a better finish due to its unique and uniform abrasive
- Does not release unwanted odors, its characteristics make it different to traditional polishers.

Advantages

- Works faster than pumice
- Contains an antiseptic that attacks the majority of bacteria & prevents the risk of cross contamination when new and old prosthetics are polished in the same equipment.



Bag with 1 lb	102-0010
Box with 25 lb	102-0011

Shinning paste for acrylics, resins, and flexible materials

polyshine

Applications

Shinning paste for acrylic denture, resins, and flexible materials to be used after Polycril

Characteristics

Paste with very abrasive agents used in low RPMs

Advantages

Small quantities are used to gain the highest permanent shine



Tube 3.5 oz	103-0020
-------------	----------

Polishing and Shinning

Rhino

Applications

Acrylics, metals and porcelain polishing paste

Characteristics

- White: All purpose
- Yellow: Metals and removes scratches (before green paste)
- Green: All purpose CrCo and CrNi – shinny finish
- Red: Gold class Au, AiCu, AlCu
- Pink: Acrylics



	Pink	White	Yellow	Green	Red
Tube 8.8 oz	104-0010	104-0050	104-0030	104-0040	104-0020

Articulator



Applications

Used for crowns and partial dentures, implants and full prosthetic models

Characteristics

- Removable magnetic mounting plate
- Incisal guide pin allows model to rest on base plate upon opening
- Colors: Blue, Pink, White, Red
- Color Chrome

Advantages

- Clear view of posterior-anterior angles
- Heavy duty
- Light weight frame
- Ergonomic design



Box w/ 1 articulator	WHITE: 123-0012	BLUE: 123-0010
Bag w/ 1 platinum plate	123-0015	
Bag w/ 10 fasteners	123-0014	

* Base plate sold separately

Economy Chrome Articulator

Applications

- Used to mount lab models

Characteristics

- Chrome finishing
- Occlusal articulator hinge that allows opening and closing

Articulator	205-0010
-------------	----------



Cyan

Impression material siliconized
with antiseptic properties



Applications

- Impressions for removable prosthesis, full skeletal and prosthodontics
- Impressions of soft tissues of the oral cavity and surrounding areas

Characteristics

- Extra fine powder
- Cherry scent
- Dimensional stability
- Dustless
- Packaged in a trilaminated metallic bag to ensure stability of the material

Advantages

- Contains chlorhexidine, an antiseptic which helps eliminate cross contamination
- Thixotropic, does not run and has a thick consistency
- Smooth texture
- High detail reproduction
- Excellent strength, flexibility and elasticity, does not tear when removed
- Detail reproduction 0.025 mm



Complies with ISO 1563/ ADA No. 18. Non-Toxic



Note: the temperature and the chemistry of the mixed water (that which is not distilled or deionized) may affect setting and working time.

Once package has been opened, seal properly to protect alginate from humidity and contamination.

	Type I	Type II
Bag of 1 lb.	101-0023	101-0026

Chromatic

Color changing impression material.



Applications

- Fixed and removable prosthesis

Characteristics

- Three phases: violet / pink / white
- Mint scent
- Extra fine powder
- Dustless Alginate
- Dimensional stability
- Packaged in a trilaminated metallic bag to ensure stability of the material

Advantages

- 130 hours without alterations or losing any details
- Good resistance
- Adequate working time
- Precision of color phase vs time
- Optimum physical qualities, fluency and consistency
- Excellent compatibility with plaster, does not interfere with setting time



Complies with ISO 1563/ ADA No. 18. Non-Toxic



Note: the temperature and the chemistry of the mixed water (that which is not distilled or deionized) may affect setting and working time.

Once package has been opened, seal properly to protect alginate from humidity and contamination.

Type I	
Bag of 1 lb.	101-0046

Orthodontic

Impression material for braces without wire



Applications

- Orthodontic use

Characteristics

- Extra fine powder
- Dimensional stability
- Banana scent
- Yellow color
- Packaged in a trilaminated metallic bag to ensure stability of the material

Advantages

- Helps avoid tearing when taking dental impressions with braces



Complies with ISO 1563/ ADA No. 18. Non-Toxic



Note: the temperature and the chemistry of the mixed water (that which is not distilled or deionized) may affect setting and working time.

Once package has been opened, seal properly to protect alginate from humidity and contamination.

Type I

Bag of 1 lb.	101-0036
--------------	----------

All Purpose

Versatile impression multipurpose material



Applications

- Fixed removable prosthesis, and full dentures

Characteristics

- Dimensional stability
- Extra fine powder
- Cherry scent
- Packaged in a trilaminate metallic bag to ensure stability of the material



Note: the temperature and the chemistry of the mixed water (that which is not distilled or deionized) may affect setting and working time.

Once package has been opened, seal properly to protect alginate from humidity and contamination.

	Type I	Type II
Bag of 1 lb.	101-0073	101-0076

Creamy impression material that is easy to use

Applications

- Removable prosthetics and skeletal

Characteristics

- White color
- Tropical scent

Advantages

- Packaged in a trilaminated metallic bag to ensure stability of the material

Note: the temperature and the chemistry of the mixed water (that which is not distilled or deionized) may affect setting and working time.

Once package has been opened, seal properly to protect alginate from humidity and contamination.

BioJel



Complies with ISO 1563/ ADA No. 18. Non-Toxic



	Type I	Type II
Bag of 1 lb.	101-0015	101-0076

Mixing Bowl

Ideal size for mixing. Flexible material that allows it work with plaster.

Blue	271-0001
Green	271-0002
Turquoise	271-0009
Black	271-0003

Sale per pack of 10 pieces



Mixing spatula

Firm yet flexible!

Blue	271-0011
Green	271-0012
Turquoise	271-0015
Black	271-0013

Individual sale



Professional Latex Dental Dam

Latex Sheet for operating area

nic tone

Isolates from one single tooth to entire dental arch.
Allows ample vision due to its contrasting color.

Helps avoid cross contamination by preventing dental material and instruments from falling into the oral cavity.

THE DENTAL
ADVISOR
++++

Characteristics

- Blue, Green or Black
- Sizes: 5" X 5" and 6" X 6"
- Three thicknesses

Advantages

- Stretches up to 800% without tearing
- Low powder usage
- Colors contrast with oral tissues
- Adjusts to the retraction of soft tissues
- Made from refined latex



	Thin	Medium	Thick
Box w/36 pcs			
Green 5" x 5"	144-0011	144-0021	144-0031
Green 6" x 6"	144-0013	144-0023	144-0033
Blue 5" x 5"	144-0010	144-0020	144-0030
Blue 6" x 6"	144-0012	144-0022	144-0032
Black 6" X 6"	-	144-0038	144-0037

	Thin	Medium	Thick
Box w/52 pcs			
Green 5" x 5"	144-0015	144-0025	144-0036
Blue 5" x 5"	144-0014	144-0024	144-0034

Dental Enamel Cleaning Kit

Antivet® Kit

THE DENTAL
ADVISOR
++++

It acts directly on stains caused by fluorosis, smoking or drinking without damaging the enamel.

An effective and reliable dental stain remover with only one dental visit.

A Fast and permanent system.

ANTIVET is specially formulated for the cleaning of teeth with pigmentation (brownish stains) on the surface of the enamel because of dental fluorosis.

The formula based acid is chemically stabilized to a controlled pH. This is necessary as it reacts with the positive ions in the forming minerals, which clean the structure of the enamel of the teeth without affecting it. It only extracts organic and inorganic substance ions that are occupying the spaces that exist between calcified crystals of fluoride.

Advantages

- Do not harm the dental enamel
- Permanent cleaning (due to fluorosis)
- Results in the first session with the dentist
- No post occurs
- Product controlled Mexican Health Ministry (Cofepris under Reg. 1522C2013 SSA.)



Clinical cases available at
www.antivet.com


Kit	146-0010
<ul style="list-style-type: none"> - Cleaning Solution 0.34 oz. - Neutralizing Solution 0.34 oz. - 2 brushes - 1 godete - 5 dental dam sheets - 1 Instructions 	



FLUOROSIS STAIN



SMOKING STAIN



Powder adhesive for dental prosthesis

- Jar with flip top lid makes application easier and prevents openings being clogged due to moisture
- Easy to apply
- Prevents the proliferation of bacteria causing bad breath
- Thanks to its hydro-active formulation it produces a protective seal that prevents food accumulation inside the prosthesis
- It is easy to clean without leaving residues
- Adhesion of up to 12 hours



Jar 0.7 oz	147-0012
------------	----------

Adhesive paste for dental prosthesis

- Up to 12 hour hold with only one application
- It allows you to chew food naturally
- It avoids bacteria proliferation that cause bad breath
- Due to its hydro-active formulation it produces a protective seal preventing food remains from accumulating inside the prosthesis
- It makes for easy cleaning without leaving food residues
- Stability for up to 12 hours



Tube 1.4 oz	147-2010
Tube 2.6 oz	147-2030

Cleansing effervescent tablets

- For dental prosthesis and ortodontic appliance. Always clean and with a nice taste!
- It eliminates 99.9% of bacteria that cause bad breath
- It keeps the prosthesis stain free
- No bad odor
- It cleans the most inaccessible areas
- Confident dental container with a small drain basket and seal system is recommended



Tube with 20 pills	147-0020
--------------------	----------

Denture container

- Discrete and practical design in order to maintain dental prosthesis
- It protects against pets
- Very practical during trips
- Easy to use
- Manufactured with antibacterial plastic
- It includes a small basket for rinsing which makes cleaning easier



Denture container	147-0030
-------------------	----------

Denture brush

- Works on dentures, partials, retainers & sport protectors
- Reach unaccessible areas
- Removes stains
- Ergonomic design
- Lasts up to 3 times longer than a regular dental brush



Denture brush	147-3010
---------------	----------

Inflex

Automatic Injection machine
for thermoplastic materials

Get a machine that will give you consistent
and accurate results.


Inflex is used to melt TCS thermoplastic material
and inject it into the flask with an easy operation,
long lifespan, low maintenance, the security of
success every time.

Characteristics:

- Digital Control
- 3 Different programs
- Pneumatic Pressure until 10 k of air.
- Stainless steel
- Adjustable Height
- 250w Resistance
- 120v

Advantages:

- LCD Screen indicate temperature, time,
type of program and actual phase of the process.
- Temperature reaches up to 662° F (350° C)
- Fast heating
- Optimum injection process
- Prevents material color changes with consistent temperature
- Equal Pressure





USA

17800 South Main Street, Suite 201
Gardena CA, 90248
Office: (310) 352 6980
Fax: (310) 352 6982
Toll free: (877) 632 6377
info@mdcdental.com

Guadalajara-Mexico (*Internacional Headquarters*)

Industria del Plástico 2113
Fracc. Industrial Zapopan Norte
Zapopan, Jalisco, México, C.P. 45130
01 800 3637 800
mdc@dental.com.mx

Certified for international quality products



Reg. ES-GS-0015/2016



Reg. ES-0472/2016



www.mdcdental.com