

+ innovo™ puraclav



+ innovo™ puraclav

COMPACT

B CLASS AUTOCLAVE

THE AUTOCLAVE OF CHOICE FOR BEAUTY, NAILS, TATTOO
AND STERILISATION OF SMALL INSTRUMENTS



AN EASY TO USE & COST EFFECTIVE STERILISATION SOLUTION


- ✓ Safe
- ✓ Compliant
- ✓ Reliable
- ✓ Easy To Use
- ✓ Cost Effective

The Innovo Puraclav Compact B Class Autoclave ticks all the boxes.

The Australian Designed Innovo Puraclav Compact is designed to be a cost effective solution that is TGA approved and compliant to AS/NZS4815 & AS/NZS4187.

Features & Benefits

- Best return on investment (low initial cost, guaranteed reliability & service)
- Low cost servicing and validations
- Robust and quality build
- B Class Autoclave - highest level of sterilisation making it suitable for use on hollow instruments and multiple wrapping
- 4 Preset Sterilisation Cycles included operating from 121deg to 134deg C including solid, wrapped and hollow instruments cycles
- LED Display with preset buttons
- B&D/Helix/Vaccum Test cycle presets



What's included

- Innovo Compact Autoclave
- Built in Thermal Printer
- 3 trays
- 1 year Parts & Labour warranty

Optional Extras

- Installation & Initial Validation
- Onsite Training
- Extended Warranty - extend the warranty out to 3 years for all parts and labour
- 4L Water Distiller - Autoclaves require the use of distilled water. The 4L Water Distiller turns regular water into distilled water ready to use in your autoclave

Specifications

Chamber Diameter:	Tray Size:
8L - 17cm	8L - 12.5cm x 25.5cm
12L - 20cm	12L - 15cm x 28cm
Chamber Depth:	Net Weight:
8L - 32cm	8L - 45kg
12L - 36cm	12L - 47.7kg

Screen Type - LED
Overall Size - 56cm x 44cm x 36cm
Designed to suit standard 500mm benchtop

AUSTRALIAN DESIGNED,
TGA APPROVED,
COMPLIANT TO WORKSAFE
& AUSTRALIAN STANDARDS



INTEREST FREE FINANCE
OPTIONS AVAILABLE





Ph 1300 651 958

Web www.nesgroup.com.au

Email enquiries@nesgroup.com.au