



FM-200

Doppler Fetal Heart Rate Detector



FM-200 Doppler Fetal Heart Rate Detector



Water-proof Probe



Back light LCD



Operation buttons



Features:

- Handheld & elegant design;
- Back light LCD;
- Display of FHR, storage, work model, battery status, signal indication and probe frequency;
- Three work modes: real time, average and manual;
- FHR alarm limits is adjustable and can be saved;
- Low power consumption design;
- Unique heart rate algorithm enable FHR more accurate and reliable;
- Buckle design for easy take;
- Automatically power off while no signal within 2 minutes;
- Inter-changeable water-proof probe of 2MHz or 3MHz.

Options:

3MHz probe, Earphone, rechargeable battery.

Specifications:

Ultrasonic Work Frequency: 2MHz, 3MHz(Optional);
Heart Rate Range: 50bpm~240bpm;
Alarm Limit: Adjustable lower limit range: 50bpm~120 bpm;
Adjustable upper limit range:160bpm~240 bpm;

Power Requirement: 2 x 1.5V (AA) Alkaline Battery or Rechargeable Battery;
Dimension: 130mm×89mm×34mm (L×W×H);
Net Weight: 405g (including battery).