

+ innovo™ puraclav

+ innovo™ puraclav

PREMIUM

B CLASS AUTOCLAVE

DESIGNED IN AUSTRALIA FOR QUALITY & RELIABILITY



AN EASY TO USE & COST EFFECTIVE STERILISATION SOLUTION

- ✓ Safe
- ✓ Compliant
- ✓ Reliable
- ✓ Easy To Use
- ✓ Cost Effective

The Innovo Puraclav B Class Autoclave ticks all the boxes.

The Australian Designed Innovo Puraclav is designed to be a reliable and cost effective solution that is TGA approved and compliant to AS/NZS4815 & AS/NZS4187.



Features & Benefits

- Best return on investment (low initial cost, guaranteed reliability & service)
- Robust and quality build
- High Quality German brand Thomas Vacuum Pumps & ULKA Water Pumps
- B Class Autoclave - highest level of sterilisation making it suitable for use on hollow instruments and multiple wrapping
- 10 Preset Sterilisation Cycles included operating from 121deg to 134deg C including solid, wrapped and hollow instruments cycles
- LED Display with preset buttons
- B&D/Helix/Vacuum Test cycle presets
- Standard with both built In Printer & USB



What's included

- Innovo Puraclav Premium Autoclave
- Built in Thermal Printer & USB
- 3 trays
- 2 year Parts & Labour warranty

Optional Extras

- Initial Validation
- Onsite Training & Installation
- Extended Warranty - extend the warranty out to 3 years for all parts and labour
- 4L Water Distiller - Autoclaves require the use of distilled water. The 4L Water Distiller turns regular water into distilled water ready to use in your autoclave


Specifications

Chamber Diameter:	Chamber Depth:
12L - 20cm	12L - 36cm
18L - 25cm	18L - 36cm
23L - 25cm	23L - 45cm

Tray Size:	Net Weight:
12L - 15cm x 28cm	12L - 47.7kg
18L - 19cm x 28cm	18L - 50.25kg
23L - 19cm x 38cm	23L - 53.5kg

Overall Size:
63.1d x 43.6w x 48.8h - (Position of feet designed to suit standard 500mm benchtop)

AUSTRALIAN DESIGNED,
TGA APPROVED,
COMPLIANT TO WORKSAFE
& AUSTRALIAN
STANDARDS



THOMAS



Recommended with Purchase



Helix Test System

- This test is suitable for every sterilisation process in a loaded autoclave. It can also be used to ensure maximum steam penetration where hollow instruments, or instruments with small Lumina are concerned.

Medicom SafeSeal Quattro Pouches

- SafeSeal Quattro pouches give the user the assurance of confirming sterilisation parameters in one easy step: time, temperature & steam with internal & external indicators printed directly on the pouch. All Medicom SafeSeal pouches provide peace of mind, convenience and economy. Properly stored, they maintain sterility of the instruments for up to 1 year post sterilisation.



Automatic Water Distiller

- Produces 4L of distilled water in 7 hours
- Easy to clean & maintain
- Automatically shuts off when finished

AUSTRALIAN DESIGNED, TGA APPROVED, COMPLIANT TO WORKSAFE & AUSTRALIAN STANDARDS

**INNOVO PURACLAV
PREMIUM AUTOCLAVE**

nesgroup

ne