

iH-B ABP Recorder is lightweight, convenient to use, and energy-efficient, making it ideal for everyday use, the recorder employs an advanced digital signal processor and a flash memory to be able to record up to 500 BP measurements.



# **iH-B** Ambulatory Blood Pressure Recorder

# iH-B

## Ambulatory Blood Pressure Recorder

### Features (Recorder):


- Check BP result on the LCD display;
- Pressure protection handled by double sensors to meet the requirement of IEC 60601-2-30;
- Storage medium: Flash memory;
- Light weight, and energy-efficient.




Easy to use, reliable, fast and well-arranged – iH-B analysis software allows transfer BP measurement data from the Recorder to a Windows, PC-based computer program for the purpose of creating reports and printouts. Physicians can review artery pulsation waveform.

### Features(Software):


- Quick and accurate analyzing system;
- Reviewable analysis techniques to obtain BP measurement ;
- Highly intellectualized airflow control;
- introduction of advanced AGC (auto-gain-control) technology and components with a large
- dynamic range, which automatically adjusts to the variation of the patients' pulsation signal range.



Blood pressure histograms



Blood pressure & pulse rate trends



Reviewable artery pulsation waveform & tabulated numerical data results

Distributed by:



PO Box 530, Camden NSW 2570

Ph: 1300 651 958

Email: enquiries@nesgroup.com.au

www.nesgroup.com.au