



B. BRAUN ADHESIVE REMOVER

SILICONE BASED ADHESIVE REMOVER

B. Braun Wound Care. Right. From the start.

B. BRAUN ADHESIVE REMOVER



HOW TO APPLY:

Spray the adhesive remover around the pouch and wait a few moments to let the adhesive work. Continue to use it as you remove your pouch. If it is not easy to remove, do not pull. Add more adhesive and wait a little before trying to pull it again. The spray dries very quickly, once dried it will not affect further appliances.



SPECIAL NOTE:

Advantages thanks to the Bag-on-valve technology:

Safe & pure: material & propellant are separated

Practical: can be used with the canister in any position and stored horizontally without leakage

Environmental benefits: use an eco-friendly propellant, nitrogen and fully recyclable.

Discreet: light noise of the spray

Longer shelf life : 3 years, substances contained by a four layer vacuum bag

Economic: no waste, can be use until it is completely empty

*data file

ADVANTAGES*:

- Does not sting on application
- Alcohol and odour free
- Preserve skin integrity
- Help the removal of the adhesive and minimise the negative impact on the skin
- Does not affect adhesion of future dressings
- Silicone based spray

COMPOSITION:

100% Siloxanes (HMDSO)

INTENDED USE:

Helps to remove Ostomy devices without pain.

For fast and painless removal of adhesives and residues from the skin surface.



PRODUCT INFORMATION:

Product Description	Size	Product Code
Adhesive remover spray	50 ml	1005-1