

# LABCON CENTRIFUGE TUBES



Labcon has a centrifuge tube for every application!



#### TUBE CATEGORIES

Metal-free, UV-blocking, thermochromic, high-speed centrifugation, and standard centrifugation.



#### VALIDATED QUALITY

Products are lot tested to ensure quality, purity, and performance. Registered to ISO 9001:2015 quality management system.



#### SUSTAINABLE PACKAGING

Designed for a low environmental impact, incorporating renewable and/or recyclable materials. Plastic racks can be autoclaved or recycled.

With the highest certified spin rates in the industry and extensive quality control testing, you can trust Labcon centrifuge tubes to reliably meet the needs of your experiment. For safe centrifugation, the Labcon quality control team performs stress tests beyond product specifications to ensure that our centrifuge tubes will meet the stated specifications in your experiments. Labcon also tests for volume line accuracy, tube wall thickness, concentricity, clarity, and leak-resistance to challenging chemicals. Additional laboratory testing confirms all products are free of biological contamination to ensure repeatability in your experiments. Labcon centrifuge tubes are made of high-quality, lab-certified, non-cytotoxic, 100% medical-grade resin and printed using low volatility, UV-fixed inks that resist organic solvents. Tubes are manufactured, packaged, and sterilized in the United States, with many tube offerings achieving ACT score status for reduced environmental impact. Labcon is a US FDA registered medical device facility and offers sterile products validated for medical, clinical and research use.

- **Wide selection of centrifuge tube sizes**

*0.5mL, 1.5mL, 1.7mL, 2.0mL, 5mL, 15mL, and 50mL*

- **Centrifuge tubes are certified free of:**

*RNase/DNase, ATP, Endotoxin/Pyrogens (<0.05EU/mL), Cytotoxic effects, and trace heavy metals (<1ppm)*

- **Sterile tubes are validated ISO 11137 for research, clinical and diagnostic use**

*Certified sterile with a 10<sup>-6</sup> sterility assurance level (SAL)*

- **Leak-resistant flat caps with sealing rings**









*Available exclusively with SuperClear® and ViewPoint™ tubes, these leak-resistant flat caps feature an elastomeric sealing ring molded securely in the cap. This ring forms a liquid-tight seal and protects from leakage of dangerous chemicals.*

- **Quality certification and ACT scores available**

*Certificates outlining functional and biological testing are available for each lot number. Specific products include an ACT environmental impact score.*

# Labcon Centrifuge Tubes

 This icon indicates products independently audited for the ACT environmental impact responsibility score. See the product score at: [act.mygreenlab.org/consumables.html](http://act.mygreenlab.org/consumables.html) 

Part Number	Product Description	RCF	Cap	Sterile	ACT	Packaging	Pack	Case
<b>SuperClear®, SuperSlik®, and ViewPoint™ Microcentrifuge Tubes</b>								
3022-800-000-9	0.5mL SuperClear® with boil-proof caps	35,000	Snap			Resealable Bags	1,000	10,000
3025-570-000-9	0.5mL SuperSlik® low protein binding surface	35,000	Snap			Resealable Bags	500	5,000
3016-870-000-9	1.5mL SuperSpin™	40,000	Snap			Resealable Bags	500	5,000
3016-871-000-9	1.5mL SuperSpin™ amber for light sensitive samples	40,000	Snap			Resealable Bags	500	5,000
3039-560-000-9	1.5mL SuperSlik® low protein binding surface	20,000	Snap			Resealable Bags	250	2,500
3005-840-000-9	1.5mL ViewPoint™ thermochromic tube	40,000	Snap			Resealable Bags	500	1,000
3034-335-000-9	1.7mL SuperClear®	15,000	Snap	✓		Resealable Bags	500	5,000
3012-870-000-9	1.5mL SuperClear® with boil-proof caps	20,000	Snap			Resealable Bags	500	5,000
3046-870-000-9	2.0mL SuperClear® with boil-proof caps	20,000	Snap			Resealable Bags	500	5,000
3046-871-000-9	2.0mL SuperClear® amber for light sensitive samples	20,000	Snap			Resealable Bags	500	5,000
3048-860-000-9	5mL SuperClear®	25,000	Snap			Resealable Bags	250	2,500
3048-861-000-9	5mL SuperClear® amber for light sensitive samples	25,000	Snap			Resealable Bags	250	2,500
<b>SuperClear® 15 and 50mL Centrifuge Tubes for High Speed Centrifugation</b>								
3131-345-008-9	15mL SuperClear® with sealing ring in cap	17,000	Flat	✓		Resealable Bags	50	500
3141-335-018-9	15mL SuperClear® with sealing ring in cap	17,000	Flat	✓		Cardboard Rack	50	500
3141-335-028-9	15mL SuperClear® with sealing ring in cap	17,000	Flat	✓		Plastic Rack	50	500
3191-335-018-9	50mL SuperClear® with sealing ring in cap	20,000	Flat	✓		Cardboard Rack	50	500
3191-335-028-9	50mL SuperClear® with sealing ring in cap	20,000	Flat	✓		Plastic Rack	25	500
3181-345-008-9	50mL SuperClear® with sealing ring in cap	20,000	Flat	✓		Resealable Bags	50	500
<b>PerformR® 15 and 50mL Centrifuge Tubes For Standard Centrifugation</b>								
3136-345-008-9	15mL PerformR® with standard cap	12,500	Flat	✓		Resealable Bags	50	500
3146-335-018-9	15mL PerformR® with standard cap	12,500	Flat	✓		Plastic Rack	50	300
3186-345-008-9	50mL PerformR® with standard cap	12,500	Flat	✓		Resealable Bags	50	500
3196-335-018-9	50mL PerformR® with standard cap	12,500	Flat	✓		Plastic Rack	25	300
<b>PerformR® 50mL Freestanding Centrifuge Tubes</b>								
3095-340-008-9	50mL PerformR® with standard cap	12,500	Flat			Resealable Bags	50	500
3095-345-008-9	50mL PerformR® with standard cap	12,500	Flat	✓		Resealable Bags	50	500
<b>UVsafe™ 15 and 50mL Safe Amber Centrifuge Tubes for Light-Sensitive Samples</b>								
3132-445-008-9	15mL UVsafe™ amber for light sensitive samples	12,500	Plug	✓		Resealable Bags	50	500
3142-435-008-9	15mL UVsafe™ amber for light sensitive samples	12,500	Plug	✓		Cardboard Rack	50	500
3182-445-008-9	50mL UVsafe™ amber for light sensitive samples	12,500	Plug	✓		Resealable Bags	50	500
3192-435-008-9	50mL UVsafe™ amber for light sensitive samples	12,500	Plug	✓		Cardboard Rack	25	500
<b>Metal Free™ 15 and 50mL Centrifuge Tubes For Applications Requiring The Assurance Of Low Background Metals</b>								
3134-345-001-9	15mL Metal Free™ with standard cap	17,000	Flat	✓		Resealable Bags	50	500
3144-335-001-9	15mL Metal Free™ with standard cap	17,000	Flat	✓		Cardboard Rack	50	500
3184-345-001-9	50mL Metal Free™ with standard cap	20,000	Flat	✓		Resealable Bags	50	500
3194-335-001-9	50mL Metal Free™ with standard cap	20,000	Flat	✓		Cardboard Rack	25	500
<b>ViewPoint™ 15 and 50mL Thermochromic Centrifuge Tubes</b>								
3210-320-008-9	15mL ViewPoint™ with sealing ring in cap	12,500	Flat			IntegraPack®	50	100
3210-335-008-9	15mL ViewPoint™ with sealing ring in cap	12,500	Flat	✓		Plastic Rack	50	300
3220-320-008-9	50mL ViewPoint™ with sealing ring in cap	12,500	Flat			IntegraPack®	50	100
3220-335-008-9	50mL ViewPoint™ with sealing ring in cap	12,500	Flat	✓		Plastic Rack	25	300

