



White Finish Double Buret Clamp

The VEE GEE white finish double buret clamp is made of aluminum, die casted and painted white.

Black Rubber sleeves allow for firm gripping of burets. The clamps are adjustable and hold two burets from micro to 100mL sizes.

The thumb style screw key attaches easily to support rods.



Specifications

Cat. No.	Capacity (mL)	Style	Pack Qty
3021	For: Micro to 100mL	Double, White Finish	1 (Each)



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com

No. VGBRC030116-3021
Printed in the USA