

VEE GEE®

VEE GEE Model 3040 Alcohol Burners feature a sturdy, heavy-walled glass design with a large capacity of 150ml. Fitted with a nickel-plated brass cap used to extinguish the flame and to prevent the evaporation of alcohol. Either methyl or ethyl alcohol can be used with the lamp.

Wick:
Long-burning, easily replaceable.

Wick Holder:
Durable, nickel-plated brass.

Burner:
Large 150ml capacity;
heavy-walled glass.

Cap:
Durable, nickel-plated brass.

- Glass Lamp w/ Nickel-Plated Brass Cap
- Large 150ml Capacity
- Ideal For Students

Alcohol Burners Specifications

Cat. No.	Description	Capacity	Diameter	Height	Weight
3040	Glass Lamp w/Nickel-Plated Brass Cap	150ml	80mm	100mm	240g
3040-W	Replacement Wicks (Pk/10)				



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com

No.VGBRC082213/3040
Printed in the USA