

J A A R



Hydrometer Jar with 1000mL ASTM Calibration Mark

The hydrometer jar used in Testing to ASTM D422 standards requires a calibrated graduation mark at 1000mL. Our VEE GEE brand hydrometer jar, of borosilicate glass is widely used for Soil Testing and other applications requiring the ASTM specifications set for hydrometer readings.

No other markings on the jar allow for a clear reading of the hydrometer and temperature while submerged within the sample liquid. The permanent enameled markings are blue and easy to see against the internal sample content.

Specifications

Cat. No.	Cat. No. Kit Components
21000-1000-L	1000mL Hydrometer Jar with ASTM Graduation Line



VEE GEE Scientific LLC

620 Lakeview Parkway • Vernon Hills, IL 60061
(800) 423-8842 • www.veege.com • sales@veege.com

VGBRHydrometerJar21000-1000-L
Printed in the USA