



Serological, Wide-Tip Pipets

Reusable Glass, For Measuring


SIBATA wide-tip-series serological pipets for biochemical and other applications. Markings and color coding are permanently enameled to the glass. SIBATA pipets are meticulously manufactured for tolerances and accuracy. Compared to plastic, reusable glass pipets are autoclavable and sustainable with less waste. Glass is a superior material for measuring samples with a wide variety of chemical compounds.

Three styles available:

- Wide-Tip, Two Lengths: 205mm-350mm
- Large Hole Opening, Short, 205mm Length
- Graduated to the Tip, 435-366mm Lengths

Specifications

Cat. No.	Capacity (mL)	Graduation Interval (mL)	Hole Diameter x Length (mm)	Color Band	Pack Qty.	Case Qty.
WideTip, Two Lengths						
20200-11A	1	0.2	1.5 x 350	Red	10	90
20200-31A	3	0.5	2.0 x 350	Green	10	90
20200-51A	5	0.5	2.0 x 350	Blue	10	60
20200-101A	10	1	3.0 x 350	Orange	10	120
20200-12A	1	0.2	1.5 x 205	Red	10	90
20200-32A	3	0.5	2.0 x 205	Green	10	90
20200-52A	5	0.5	2.0 x 205	Blue	10	60
20200-102A	10	1	3.0 x 205	Orange	10	120
Large Hole, Short Length						
20210-52A	5	0.5	3.0 X 205	Blue	10	60
20210-102A	10	0.2	3.5 X 205	Orange	10	120
Graduated to the Tip						
20140-1A	1	0.01	1.2 X 360	N/A	10	60
20140-2A	2	0.02	1.2 X 360	N/A	10	60
20140-5A	5	0.05	1.2 X 366	N/A	10	60
20140-10A	10	0.1	1.4 X 366	N/A	10	60
20140-20A	20	0.1	1.5 X 435	N/A	10	60



WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



20140-10A
1.4mm Hole Opening,
Graduated to Tip,
366mm Long



20200-101A
Wide Tip: 3.0mm Hole Opening,
350mm Long



20210-102A
Extra Wide 3.5mm Hole Opening,
Short 205mm Long