



VEE GEE® Brand 3050 Series Support Jacks are an essential companion to any laboratory application requiring lifting and lowering of instruments or apparatus. Ideal for supporting glassware, stirrers, baths, hotplates, etc. Constructed of sturdy anodized aluminum with stainless-steel crossbars; platform height is quickly and smoothly adjusted by oversized knobs. An optional 8 x 8" plate set can be attached to the jack to provide additional support and for accommodating larger instruments.

- Exceptionally Durable Construction
- Oversized Knobs For Simple Lifting
- Optional Support Plates & Rods



Support Jacks

VEE GEE®, Anodized Aluminum, Vertical Lift 2.75" to 10.5"

Cat. No.	Description	Pack Qty.
3050	Support Jack, 140 x 115 mm (5.5 x 4.5 in)	1 each
3050-200	Support Plate Set, 200 x 200 mm (8 x 8 in)	1 set of 2
3050-R305	Support Rod, Attachment, 305 mm (12 in)	1 each



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com

No.VGBRC081913/3050
Printed in the USA