

Standard Surgical Gown, 45gsm

Sterile, Single-Use

REF: BKM-GS011

UDI: 5060586271879

Size: Large (L)

Contents:

- 1 x High Performance Surgical Gown
- 2 x Hand Towels
- 1 x Double SMS Wrap



EN13795-1:2019

Instructions for Use:

Place pack on surface; peel off bag and remove gown; hold from top and shake down; place hands into sleeves; seek assistance to draw back from chest and fasten neck and waist belts and secure with tape.

CE 0197



OCT 2020



OCT 2023

LATEX FREE

LOT: 123651

DO NOT USE IF PACK IS DAMAGED OR OPEN

Ⓢ SINGLE-USE

EC REP

GENIMPEX S.R.L.S
Viale Lombardia 143, 20816,
Brugherio, Italy

Manufactured for BlueKit Medical Ltd by:
Hefei C&P Nonwoven Products Co.,Ltd
No. 22 Park Road, Feidong new city development
area, Hefei, Anhui, China

BlueKit Medical Limited, Surrey, England.
Company Registration No. 08097913

Standard Surgical Gown, 45gsm

Sterile, Single-Use

REF: BKM-GS010

UDI: 5060586271862

Size: Medium (M)

Contents:

- 1 x High Performance Surgical Gown
- 2 x Hand Towels
- 1 x Double SMS Wrap



EN13795-1:2019

Instructions for Use:

Place pack on surface; peel off bag and remove gown; hold from top and shake down; place hands into sleeves; seek assistance to draw back from chest and fasten neck and waist belts and secure with tape.

CE 0197



OCT 2020



OCT 2023

LATEX FREE

LOT:

123651

DO NOT USE IF PACK IS DAMAGED OR OPEN



SINGLE-USE

EC REP

GENIMPEX S.R.L.S
Viale Lombardia 143, 20816,
Brugherio, Italy

Manufactured for BlueKit Medical Ltd by:
Hefei C&P Nonwoven Products Co.,Ltd
No. 22 Park Road, Feidong new city development
area, Hefei, Anhui, China

BlueKit Medical Limited, Surrey, England.
Company Registration No. 08097913