




PPE Treatment Kit

Product Specification Sheet

Isolation Gown

Single-Use, Isolation-Gown with Knitted Cuff

Material	PP with PE coating
Weight	30gsm: 18gsm PP + 12gsm PE +2g
Size	M size: 120 x 140 cm L size: 130x150cm
Colour	Blue
Standard	CE EN14126
Application	General Medical, Labs, Patient Care, and Critical Care Units



Neck tie, around back



Long sleeve, knitted cuffs



Tie around the waist

Our all-in-one PPE Kit provides all the protection you need to practice safely and efficiently.

Reduce risk, instill confidence and increase efficiency.

Please note that VAT is temporarily zero-rated on the sale of PPE items until the 31st October.

T: +44 (0) 20 8224 2883 E: info@bluekitmedical.com W: www.bluekitmedical.com
BlueKit Medical | Capital Business Centre | South Croydon | CR2 0BS | United Kingdom
Company No: 08097913 VAT No: 141 0046 73

PPE Treatment Kit
Product Specification Sheet


Isolation Gown
Single-Use, Isolation-Gown with Knitted Cuff

Medium (M)
Size Specification:

A	B	C	D	E	F
120cm +3cm	140cm +3cm	22cm +1cm	57cm +2cm	28cm +1cm	16cm +1cm
G	H	I	J	K	M
35cm +3cm	95cm +3cm	7cm +1cm	7cm +1cm	5.5cm +0.5cm	62cm +3cm

Large (L)
Size Specification:

A	B	C	D	E	F
130cm +3cm	150cm +3cm	22cm +1cm	57cm +2cm	28cm +1cm	16cm +1cm
G	H	I	J	K	M
35cm +3cm	95cm +3cm	7cm +1cm	7cm +1cm	5.5cm +0.5cm	62cm +3cm



Our all-in-one PPE Kit provides all the protection you need to practice safely and efficiently.

Reduce risk, instill confidence and increase efficiency.

Please note that VAT is temporarily zero-rated on the sale of PPE items until the 31st October.