

Medintim



**Active<sup>3</sup>**  
Erection System

**Vacuum therapy for Erectile Dysfunction (ED)**

Medintim

Active<sup>3</sup><sup>®</sup>  
Erection System

Medintim

Place the  
guide sleeves  
and the cone.

Active<sup>3</sup>  
Erection System

Medintim

Slide on a  
constriction ring  
using the cone.

Active<sup>3</sup>  
Erection System



Medintim

Alternatively,  
slide on a  
constriction ring  
using your hands.

Active<sup>3</sup>  
Erection System

Medintim

Ready for  
application.



Active<sup>3</sup>  
Erection System

A person is shown from the side, applying a blue patch to their upper arm. The person's hand is holding the patch, and another hand is visible near the application site. The background is a solid blue color.

Medintim

**Application on the body**



Medintim

Apply  
lubricant.

Active<sup>3</sup>  
Erection System




Medintim



Hold the device  
against your body  
and pump.

Active<sup>3</sup>  
Erection System

Medintim



The vacuum in the cylinder creates an erection.

Active<sup>3</sup>  
Erection System

Medintim



Push the  
constriction ring  
onto the penis.

Active<sup>3</sup>  
Erection System

Medintim

The constriction ring maintains the erection (max. 30 min.)

Active<sup>3</sup>  
Erection System



Medintim

**After application:  
Slip off the constriction ring  
and remove it.**


**Active<sup>3</sup>**  
Erection System

A photograph of a man's hands and lower torso. He is holding a clear plastic ring in his right hand and a small, dark, cylindrical object in his left hand. A blue callout box with a white border and a blue dot points from the text 'State after application.' to the man's penis, which is partially visible and appears to be in a state of erection. The background is a warm, orange-brown color.

Medintim

State after application.

Active<sup>3</sup>  
Erection System



[www.medintim.de](http://www.medintim.de)

CE 0483

FDA registered device  
Made in Germany

# Medintim

[www.medintim.de](http://www.medintim.de)



KESSEL Medintim GmbH  
Kelsterbacher Str. 28  
64546 Mörfelden-Walldorf  
Germany  
+ 49 (0) 6105/203 720