

PA-BED / PA-CHAIR



PA-BED / PA-CHAIR

Product Description

The alarm systems contain the alarms and sensor pads for bed or chair. The systems are part of Fall Prevention Program. They do not intend to replace the protocols and procedures advised by Fall Prevention Program guidelines.

Caution: This product is not a restraint. It will not prevent a resident from getting out of their bed or chair. This product is designed as a monitor, alerting the staff when the resident has risen from their bed or chair.

This product is to be used in conjunction with a total fall prevention program

Clean Instruction

- Clean the alarm pad using a surface disinfectant only
- DO NOT sterile the pad in high temperature



Specification

Feature	Tamper-Resistant Alarm	Universal Alarm
Nursing Call (optional)	V	X
Volume Adjustable	V	X
Tempo Adjustable	V	X
Tone Adjustable	V	X
	Bed	Chair
Sensor Pad Size	75.5cm x 30cm(±3%)	30cm x 30cm(±3%)

