

PARAMOUNT BED CO., LTD. 2-14-5 Higashisuna Koto-Ku, Tokyo 136-8670, Japan Tel: +81 3 3648 2961 Fax: +81 3 3648 2420 http://www.paramount.co.jp/english/

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.





one small step from the bed, a giant leap for their lives

Reduces Risks

Prevention of Patient Falls

It is important to keep a firm footing on the floor when the patient is getting up from the bed. A low height bed, allows even users with small physiques to sit on the side of the bed with their heels firmly on the floor. It also helps to reduce the impact in the event of a fall.

Prevention of Pressure Ulcers

During back rest elevation, the back of the frame is extended to accommodate the natural lengthening of the body as the patient's back is raised. This movement not only prevents the patient from sliding down the mattress toward the foot side of the bed, but also reduces the risk of pressure ulcers.



Infection Control

Flat surface design of the headboard, footboard and side-rails allows for easier cleaning.

Promotes Rehabilitation

<Supports Rehabilitation in Stages>

Step 1: Shifting of Sleeping Positions Step 2: Sitting Up with Assistance

Step 3: Training in Sitting on the Side of the Bed Step 4: Training in Standing Up

Others

Manual CPR Levers Attachment Parts for Restraint Belts Angle Indicators (2 Angle Indicators / Bed)

Additional Accessories

IV Pole Lifting Pole Oxygen Cylinder Holder (Hold up to 500L or 1,000L) Extension Frame & Mattress Bed Linen Rack Battery Traction Frame



– Cabinet (PF-1201)




PATIENT CONTROL PANEL-

Color Options

ANGLE INDICATOR -

MANUAL CPR LEVER



Specifications

Model No.	PA-6315	PA-6325	PA-6318	PA-6328	PA-6319	PA-6329	PA-6310	PA-6320
Total Length	2,200mm							
Total Width	1,056mm							
Total Height	735 to 1,115mm						760 to 1,140mm	
Mattress Base Size	1,910 x 910 (mm) or 2,000 x 910 (mm)							
Mattress Base Height*	350 to 730mm						375 to 755mm	
Safe Working Load**	245kgf (2,400N)							
Caster	125mm dia. single-wheel casters with total locking system		125mm dia. double-wheel casters with total locking system		125mm dia. single-wheel casters with diagonal locking system		150mm dia. single-wheel casters with total locking system	
Back Raise Angle	0° to 70°							
Knee Raise Angle	0° to 25°							
Side-rail Integrated Controller	0		0		0		0	
Handheld Controller	—	\bigcirc	_	0		0		\bigcirc

*The mattress base height indicates the height from the floor to the mattress base surface. **Maximum load with which the bed can be operated.

(The total load of user's weight and incidental items including the mattress and options.)

X-ray / Foot Light / Anti-static Caster are also available. Please contact Paramount Bed or our distributors for more information.

CASTER OPERATION PEDAL



HANDHELD CONTROLLER



ATTACHMENT PARTS FOR RESTRAINT BELT

Medium Oak

Pink, Orange and Grand Blue are made to order.