



# Technical Sheet Ultra Slim Ostomy Belt



- Total height 27 cm with optional pocket 17cm high
- With or without pocket option.
- Custom Tailored
- Antibacterial and absorbent fabric
- Certified for medical use according to Regulation (EU) 2017/745

## Composition

85% Polyamide 15% Elastane

## LAUNDRY



Standard Machine (30°C, 86°F)



Do Not Bleach



Do Not Tumble Dry



Do Not Iron



Do Not Dry Clean

## STORAGE

Dry Place

## MANUFACTURING COUNTRY

Spain

## Customer Service

E-mail: [info@ostocare.com](mailto:info@ostocare.com)  
International : 0034 910 601 256  
Spain : 910 601 256

