

## Walk Star

Walking aid

Order no.:

-----

## Walk Star

Walking aid

Article number: RR160=0000\_K



### Standard equipment

- Walking aid positioned behind the user's body
- Folding, height-adjustable tube construction in five sizes
- With swivelling or rigidly mounted front wheels
- Models with swivelling front wheels have a swivel lock
- Rear wheels with reverse-roll lock

☐ Provide estimate only

☐ Order

Billing address/customer number

Delivery address

Sales person

Name

.....  
Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

### Frame

RR1600x00000-307	<input type="checkbox"/>	Walk Star, size ½, sparkle red, handle height 39 – 48 cm	Add. charge
RR1600x00001-307	<input type="checkbox"/>	Walk Star, size 1, sparkle red, handle height 45 – 58 cm	Add. charge
RR1600x00002-306	<input type="checkbox"/>	Walk Star, size 2, saturn blue metallic, handle height 53 – 69 cm	Add. charge
RR1600x00003-306	<input type="checkbox"/>	Walk Star, size 3, saturn blue metallic, handle height 63 – 82 cm	Add. charge
RR1600x00004-027	<input type="checkbox"/>	Walk Star, size 4, silver metallic, handle height 83 – 102 cm	Add. charge

### Front wheels

RR1600x07000	<input type="checkbox"/>	Rigidly mounted	–
RR1600x01100	<input type="checkbox"/>	Swivelling	Add. charge

### Brake

RR1610002000	<input type="checkbox"/>	Friction brake	–
RR161000x000	<input type="checkbox"/>	Driving brake and wheel lock	Add. charge

Notice: The friction brake and the driving brake with wheel lock cannot be combined with each other. Installation of the friction brake or wheel lock excludes the lowest setting position of the Walk Star size ½ and 1.

### Accessories

RR161000x00x	<input type="checkbox"/>	Vertical handles with arm supports	Add. charge
HR24100030	<input type="checkbox"/>	Pelvic support, pair	Add. charge

### Information



RR1600x01100  
Front wheels, swivelling,  
with swivel lock



RR1610002000  
Friction brake



RR161000x000  
Driving brake and wheel  
lock



RR161000x00x  
Vertical handles with arm  
supports