

## KRUUSE Rehab Carpal Protector

Designed by ReDog of Sweden™

To be used before and after surgery or in cases of chronic arthritis i.e. Osteoarthritis, OCD, ligament injuries, arthritis, ED and other assessments. When dogs suffer injuries in the lower leg joints such as from the carpal joints, there are more often traumatic causes behind this occurrence. In connection with these complex joints being injured, ligaments are often affected in some way. The healing period is long and the dog needs a support around the joint as to not worsen the injury during this period.

### User Manual

The carpal protector is stretched around the entire leg and touch fastener band in accordance with the stability one wants to achieve with the protection


### Care

Wash by hand and air-dry.



Size Chart

Size	Cat. No	Suitable for
XXX-Small	279935	Bichon Frisé, Chihuahua, Pomeranian, Yorkshire Terrier
XX-Small	279934	Boston Terrier, Cairn Terrier, Lhasa Apso
X-Small	279844	Cocker Spaniel, Whippet, Irish Terrier, Cavalier King Charles Spaniel
Small	279845	Collie, Springer Spaniel, Border Collie, Beagle
Medium	279846	Labrador Retriever, Golden Retriever, Dalmatian, Boxer, Australian Shepherd
Large	279847	Rottweiler, Bernese Mountain Dog, German Shepherd, Rhodesian Ridgeback
X-Large	279882	Large Rottweiler, Large Bernese Mountain Dog, Mastiff, Great Dane



	Size chart KRUUSE Rehab Carpal Protectors in cm	XXXS (cm)	XXS (cm)	XS (cm)	S (cm)	M (cm)	L (cm)	XL (cm)
A	Total height, from edge to edge	6.5	7.5	8.5	9.5	11.5	13.5	15.5
B	Total width, from edge to edge	8.5	9.5	10.5	12.5	14.5	17.5	20.5