

Baseline Bubble Inclinometer

Product Information



Weight	0.09 kg
Dimensions	10.2 x 0.5 x 10.2 cm
Brand	Baseline
Item No.	W50178 [1005901]

Description

The Baseline Bubble Inclinometer accurately measures complete range of motion. Measurement techniques of the bubble inclinometer have been codified and appear in the AMA Guide to the evaluation of Permanent Impairment, third edition.

Place the Baseline Bubble Inclinometer near the joint to be measured: turn dial until scale reads 0; take joint through its range; read range traveled directly from dial.

Numerous neck and back measurements require the simultaneous use of 2 inclinometers.

https://www.3bscientific.com/ph/baseline-bubble-inclinometer.p_906_6324.html