



ENSURE EVERY CPR COMPRESSION COUNTS...

NEW BEATY CPR FEEDBACK DEVICE

PART NO. 879411

Beaty was developed to provide real time CPR feedback to help first responders and bystanders perform CPR in case they do not have any prior medical knowledge.

This device could help save lives.

✓ **BEEPS WHEN COMPRESSIONS ARE EFFECTIVE**

✓ **GREAT ADDITION TO ANY FIRST AID KIT OR AED CARRY CASE**

✓ **SMART AND EASY TO USE WITH NO BUTTONS OR DISPLAYS**

✓ **CONVENIENTLY POCKET-SIZED FOR QUICK ACCESS**


SPECIFICATIONS

- Beaty runs on a 3v lithium battery that can last for more than 5 years on standby mode or 15 hours of continuous CPR
- Beaty should be used on patients 8 years and older

CONTACT US TO ORDER TODAY!

HOW TO USE BEATY:

1




In case of CPR Call for HELP

2



Place Beaty in the middle of patient's chest and Beaty will make a sound for adequate chest compressions performed, reaching a depth of 5cm.

3



Push hard until a "beep" sound is heard. If no sound is heard, this means you need to push harder

