

Radiation Protection Aprons

Drop-off Front Lead Apron

Providing comfortable front protection during short surgeries and radiology procedures



Customised
Just For YOU



Unique 3-Year
Guarantee

Product Features Overview

- Radiation Protection Standards: IEC 61331.1-2014 & IEC 61331.3-2014, AU/NZS 4543.3-2000.
- With Standard Velcro Thyroid Collar and Pocket.
- Front Protection: one panel.
- Fully adjustable Velcro strap for back support and even weight distribution.
- Padded shoulders for more comfort.
- Soft fabric material to allow for complete flexibility.

Custom Options

- Embroidery or Logo on pocket
- Bib or Cosy Thyroid collar with Magnet or Velcro closure
- Extra pocket (e.g. dosimeter)
- Access to A-Track Inventory System
- Available Fabrics:
 - Easy-Clean Flex-DC
 - Stain Resistant Spangly
 - Standard Nylon

| Protection Materials | Lead Equivalence* | Protection | Size Options |
|----------------------|-------------------|-------------------------|--|
| Lead | 0.25 | <p>Front protection</p> | XS S M L XL XXL |
| Light Lead Plus | 0.35 | | Other sizes available on request. |
| Supreme Lead Free | 0.50 | | Please refer to size chart on back page. |


* Measurement in mm Pb.

Sizing Chart for Front Lead Aprons

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.

| Front Lead Aprons (in cm) | | | | | |
|---------------------------|-----|-----|-----|-----|-----|
| Size | XS | S | M | L | XL |
| Shoulder width | 36 | 36 | 42 | 42 | 50 |
| Front width | 50 | 55 | 60 | 75 | 75 |
| Standard length | 100 | 100 | 100 | 100 | 100 |

Shoulder width 36, 42 and 50 available on request. Other lengths available: 90, 110, 120. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to info@deutschmedical.com.au and we will assist you in choosing the most suitable size for your radiation protection apron.

| My Individual Measurements | |
|--|----|
| Distance shoulder to shoulder | cm |
| Length from (highest point of) shoulders to below knee | cm |
| Width at hip | cm |